



KJM 9250

^1H - ^{15}N NMR spectra on the AVIIIHD-800 spectrometer.

Version 2.0

Topspin 3,5 Windows 7 AVIIIHD800



© Professor Emeritus Alistair Lawrence Wilkins,
University of Waikato, New Zealand.
January 2020

^{15}N - ^1H f3 (and f2) Experiments on the AVIIIHD-800 Spectrometer

1.1 Introduction

aw coded ^{15}N - ^1H f3 (and f2) HSQC, HSQC-DIPS12 and HMBC parameter sets are set up with 2048 (2K) points acquired across a 12 ppm ^1H window centered at 4.7 ppm and a 30 ppm ^{15}N window (SW) centered at 118 ppm (O3). The ^1H and ^{15}N NMR signals of peptides and microcystins typically occur in these windows. Different SW and O3 settings may (will) be required for other nitrogen containing compounds.

1.2 Processing

HSQC and HSQC-DIPS12 experiments are phase sensitive experiments which should be phased **before** using the **abs1** and **abs2** commands.

HMBC experiments are absolute value experiment. Phasing is not required.

2.0 Experiments and Parameter Sets

The following aw coded ^{15}N - ^1H f3 HSQC, HSQC-DIPS12 and HMBC parameter sets are available on the AVIIIHD-800 MHz spectrometer.

- 2.1 ^1H - ^{15}N f3HSQC spectrum**
- 2.2 ^1H - ^{15}N f3HSQCfp spectrum with a flip back pulse**
- 2.3 ^1H - ^{15}N f3HSQC135**
- 2.4 ^1H - ^{15}N f3HSQCed**
- 2.5 ^1H - ^{15}N f3TROSY**

- 2.6 ^1H - ^{15}N f3HSQC-DIPS12 spectrum**
- 2.7 ^1H - ^{15}N f3HSQC-NOESY**

- 2.8 ^1H - ^{15}N f3HMBCET spectrum**
- 2.9 ^1H - ^{15}N f3HMBCQF spectrum**

- 2.10 ^1H - ^{15}N f2HSQC**
- 2.11 ^1H - ^{15}N f2HSQC135**

- 2.12 ^1H - ^{15}N BESTf3HSQC**
- 2.13 ^1H - ^{15}N BESTf3TROSY**

2.1 ^1H - ^{15}N HSQC spectrum

Parameter sets: **awf3hsqc** (+ **getprosol**)

Pulse programme: **awf3hsqc**

Type **eda** (enter) and review the following default parameters

SW ^1H = 12 ppm, **SW** ^{15}N = **30** ppm (or other suitable values).

TD ^1H = 2K, **TD** ^{15}N = 128-160 (your choice).

NS = multiple of 8 or 16, **DS** = 16.

Type **ased** (enter) and review other parameters used in the job.

O1P = ^1H spectral window midpoint = 4.7 ppm other value of your choice.

O3P = ^{15}N spectral window midpoint = 118 ppm other value of your choice.

D1 = repetition delay = **1.0** sec or other time of your choice.

CNST4 = 1J ^{15}N - ^1H coupling constant = **90** Hz or other value of your choice.

D24 = **D26** = $1/4J$ are auto calculated from **CNST4**

ZGOPTNS = Not used.

Check gradient settings.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2)** = 2K **SI(F1)** = 512 or 1K

WDW(F1) = **WDW(F2)** = QSINE

SSB(F2) = **SSB(F1)** = 2

xfb, **abs1** and **abs2**



Expansion of the ^1H 7-9 ppm and ^{15}N 100-126 ppm region of a ^1H - ^{15}N HSQC spectrum of a peptide which has 9 amino acid units

2.2 ^1H - ^{15}N HSQC spectrum with a flipback pulse

Parameter sets: **awf3hsqcfp** (+ **getprosol**)

Pulse programme: **awf3hsqcfp**

Type **eda** (enter) and review the following default parameters

SW ^1H = 12 ppm, **SW** ^{15}N = **30** ppm (or other suitable values).

TD ^1H = 2K, **TD** ^{15}N = 128-160 (your choice).

NS = multiple of 8 or 16, **DS** = 16.

Type **ased** (enter) and review other parameters used in the job.

O1P = ^1H spectral window midpoint = 4.7 ppm or other value of your choice.

O3P = ^{15}N spectral window midpoint = 118 ppm or other value of your choice.

D1 = repetition delay = **1.0 sec** or other time of your choice.

CNST4 = 1J ^{15}N - ^1H coupling constant = **90 Hz** or other value of your choice.

D24 = $1/8J$ for all NH_x signals is auto calculated from **CNST4**

ZGOPTNS = Not used.

Check gradient settings.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K SI(F1) = 512 or 1K**

WDW(F1) = WDW(F2) = QSINE

SSB(F2) = SSB(F1) = 2

xfb, abs1 and abs2



Expansion of the ^1H 7-9 ppm and ^{15}N 110-130 ppm region of a ^1H - ^{15}N HSQCfp spectrum of a peptide which that has 9 amino acid units.

2.3 ^1H - ^{15}N HSQC spectrum

Parameter sets: **awf3hsqc135 (+ getprosol)**

Pulse programme:

Type **eda** (enter) and review the following default parameters

SW ^1H = 12 ppm, **SW** ^{15}N = **30** ppm (or other suitable values).

TD ^1H = 2K, **TD** ^{15}N = 128-160 (your choice).

NS = multiple of 8 or 16, **DS** = 16.

Type **ased** (enter) and review other parameters used in the job.

O1P = ^1H spectral window midpoint = 4.7 ppm other value of your choice.

O3P = ^{15}N spectral window midpoint = 118 ppm other value of your choice.

D1 = repetition delay = **1.0 sec** or other time of your choice.

CNST4 = 1J ^{15}N - ^1H coupling constant = **90 Hz** or other value of your choice.

D24 = **D26** = $1/4J$ are auto calculated from **CNST4**

ZGOPTNS = Not used.

Check gradient settings.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2)** = 2K **SI(F1)** = 512 or 1K

WDW(F1) = **WDW(F2)** = QSINE

SSB(F2) = **SSB(F1)** = 2

xfb, **abs1** and **abs2**



File 8 Expansion of the ^1H 7-9 ppm and ^{15}N 100-126 ppm region of a ^1H - ^{15}N HSQC spectrum of a peptide which has 9 amino acid units

2.4 ^1H - ^{15}N HSQC spectrum

Parameter sets: **awf3hsqc** (+ **getprosol**)

Pulse programme: **awf3hsqc**

Type **eda** (enter) and review the following default parameters

SW ^1H = 12 ppm, **SW** ^{15}N = **30** ppm (or other suitable values).

TD ^1H = 2K, **TD** ^{15}N = 128-160 (your choice).

NS = multiple of 8 or 16, **DS** = 16.

Type **ased** (enter) and review other parameters used in the job.

O1P = ^1H spectral window midpoint = 4.7 ppm other value of your choice.

O3P = ^{15}N spectral window midpoint = 118 ppm other value of your choice.

D1 = repetition delay = **1.0 sec** or other time of your choice.

CNST4 = 1J ^{15}N - ^1H coupling constant = **90 Hz** or other value of your choice.

D24 = **D26** = $1/4J$ are manually set

ZGOPTNS = Not used.

Check gradient settings.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2)** = 2K **SI(F1)** = 512 or 1K

WDW(F1) = **WDW(F2)** = QSINE

SSB(F2) = **SSB(F1)** = 2

xfb, **abs1** and **abs2**



File 9 Expansion of the ^1H 7-9 ppm and ^{15}N 100-126 ppm region of a ^1H - ^{15}N HSQC spectrum of a peptide which has 9 amino acid units

2.5 ^1H - ^{15}N TROSY spectrum

Parameter sets: **troseyff3gpsi** (+ **getprosol**)

Pulse programme: **awf3troty**

Type **eda** (enter) and review the following default parameters

SW ^1H = 12 ppm, **SW** ^{15}N = 30 ppm (or other suitable values).

TD ^1H = 2K, **TD** ^{15}N = 128-160 (your choice).

NS = multiple of 8 or 16, **DS** = 16.

Type **ased** (enter) and review other parameters used in the job.

O1P = ^1H spectral window midpoint = 4.7 ppm other value of your choice.

O3P = ^{15}N spectral window midpoint = 118 ppm other value of your choice.

D1 = repetition delay = 1.0 sec or other time of your choice.

CNST4 = 1J ^{15}N - ^1H coupling constant = 90 Hz or other value of your choice.

ZGOPTNS = Not used.

Check gradient settings.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2)** = 2K **SI(F1)** = 512 or 1K

WDW(F1) = **WDW(F2)** = QSINE

SSB(F2) = **SSB(F1)** = 2

xfb, **abs1** and **abs2**



File 5 Expansion of the ^1H 7-9 ppm and ^{15}N 100-126 ppm region of a ^1H - ^{15}N HSQC spectrum of a peptide which has 9 amino acid units

2.6 ^1H - ^{15}N HSQC-DIPSI spectrum

Parameter sets: **awf3hsqdieltf3gpsic** (+ getprosol)

Pulse programme: **awf3hsqc-dipsi**

Type **eda** (enter) and review the following default parameters

SW ^1H = 12 ppm, **SW ^{15}N** = **30** ppm (or other suitable values).

TD ^1H = 2K, **TD ^{15}N** = 128-160 (your choice).

NS = multiple of 8 or 16, **DS** = 16.

Type **ased** (enter) and review other parameters used in the job.

O1P = ^1H spectral window midpoint = 4.7 ppm other value of your choice.

O3P = ^{15}N spectral window midpoint = 118 ppm other value of your choice.

D1 = repetition delay = **1.0 sec** or other time of your choice.

CNST4 = $^1J^{15}\text{N}$ - ^1H coupling constant = **90 Hz** or other value of your choice.

D24 = **D26** = $1/4J$ are auto calculated from **CNST4**

D9 = **80 msec**

ZGOPTNS = Not used.

Check gradient settings.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2)** = 2K **SI(F1)** = 512 or 1K

WDW(F1) = **WDW(F2)** = **QSINE**

SSB(F2) = **SSB(F1)** = 2

xfb, **abs1** and **abs2**



File 14 Expansion of the ^1H 7-9 ppm and ^{15}N 100-126 ppm region of a ^1H - ^{15}N HSQC - DIPSI2 spectrum of a peptide which has 9 amino acid units

2.7 ^1H - ^{15}N HSQC-NOESY spectrum

Parameter sets: **awf3hsq-noesy + getprosol)**

Pulse programme: **hsqcetf3gpno**

Type **eda** (enter) and review the following default parameters

SW ^1H = 12 ppm, **SW** ^{15}N = **30** ppm (or other suitable values).

TD ^1H = 2K, **TD** ^{15}N = 128-160 (your choice).

NS = multiple of 8 or 16, **DS** = 16.

Type **ased** (enter) and review other parameters used in the job.

O1P = ^1H spectral window midpoint = 4.7 ppm other value of your choice.

O3P = ^{15}N spectral window midpoint = 118 ppm other value of your choice.

D1 = repetition delay = **1.0 sec** or other time of your choice.

CNST4 = 1J ^{15}N - ^1H coupling constant = **90 Hz** or other value of your choice.

D24 = **D26** = $1/4J$ are auto calculated from **CNST4**

D9 = **80 msec**

ZGOPTNS = Not used.

Check gradient settings.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K SI(F1) = 512 or 1K**

WDW(F1) = WDW(F2) = QSINE

SSB(F2) = SSB(F1) = 2

xfb, abs1 and abs2



File 7 Expansion of the ^1H 7-9 ppm and ^{15}N 100-126 ppm region of a ^1H - ^{15}N HSQC - NOESY spectrum of a peptide which that has 9 amino acid units

2.8 ^1H - ^{15}N HMBCE spectrum

Parameter sets: **awf3hmbcet+ getprosol)**

Pulse programme: **hmbcetf3gpnd**

Type **eda** (enter) and review the following default parameters

SW ^1H = 12 ppm, **SW ^{15}N** = **30** ppm (or other suitable values).

TD ^1H = 2K, **TD ^{15}N** = 128-160 (your choice).

NS = multiple of 8 or 16, **DS** = 16.

Type **ased** (enter) and review other parameters used in the job.

O1P = ^1H spectral window midpoint = 4.7 ppm other value of your choice.

O3P = ^{15}N spectral window midpoint = 118 ppm other value of your choice.

D1 = repetition delay = **1.0 sec** or other time of your choice.

CNST4 = $^1\text{J } ^{15}\text{N}$ - ^1H coupling constant = **90 Hz** or other value of your choice.

CNST13 = **6 Hz**

ZGOPTNS = Not used.

Check gradient settings.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K SI(F1) = 512 or 1K**

WDW(F1) = WDW(F2) = QSINE

SSB(F2) = SSB(F1) = 2

xfb and xf2m + abs1 and abs2



File 16 Expansion of the ^1H 7-9 ppm and ^{15}N 100-126 ppm region of a ^1H - ^{15}N HMBCE spectrum of a peptide which has 9 amino acid units

2.9 ^1H - ^{15}N HMBCQF spectrum

Parameter sets: **awf3hmbcqw + getprosol)**

Pulse programme: **hmbcf3gpndqw**

Type **eda** (enter) and review the following default parameters

SW ^1H = 12 ppm, **SW** ^{15}N = **30** ppm (or other suitable values).

TD ^1H = 2K, **TD** ^{15}N = 128-160 (your choice).

NS = multiple of 8 or 16, **DS** = 16.

Type **asd** (enter) and review other parameters used in the job.

O1P = ^1H spectral window midpoint = 4.7 ppm other value of your choice.

O3P = ^{15}N spectral window midpoint = 118 ppm other value of your choice.

D1 = repetition delay = **1.0** sec or other time of your choice.

CNST4 = ^1J ^{15}N - ^1H coupling constant = **90** Hz or other value of your choice.

CNST13 = **6** Hz

ZGOPTNS = Not used.

Check gradient settings.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K** **SI(F1) = 512** or **1K**

WDW(F1) = WDW(F2) = QSINE

SSB(F2) = SSB(F1) = 2

xfb and **xf2m** + **abs1** and **abs2**



File 17 Expansion of the ^1H 7-9 ppm and ^{15}N 100-126 ppm region of a ^1H - ^{15}N HMBCQF spectrum of a peptide which has 9 amino acid units

2.10 ^1H - ^{15}N f2HSQC spectrum

Parameter sets: **awhsqcetgp** (+ **getprosol**)

Pulse programme: **awf2hsqc**

Type **eda** (enter) and review the following default parameters

SW ^1H = 12 ppm, **SW** ^{15}N = **30** ppm (or other suitable values).

TD ^1H = 2K, **TD** ^{15}N = 128-160 (your choice).

NS = multiple of 8 or 16, **DS** = 16.

Type **ased** (enter) and review other parameters used in the job.

O1P = ^1H spectral window midpoint = 4.7 ppm or other value of your choice.

O2P = ^{15}N spectral window midpoint = 118 ppm or other value of your choice.

D1 = repetition delay = **1.0 sec** or other time of your choice.

CNST4 = 1J ^{15}N - ^1H coupling constant = **90 Hz** or other value of your choice.

ZGOPTNS = Not used.

Check gradient settings.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2)** = 2K **SI(F1)** = 512 or 1K

WDW(F1) = **WDW(F2)** = QSINE

SSB(F2) = **SSB(F1)** = 2

xfb, **abs1** and **abs2**



File 101 ^1H - ^{15}N f2HSQC spectrum of a peptide which that has 9 amino acid units.

2.11 ^1H - ^{15}N f2HSQC spectrum

Parameter sets: **awf2hsqc135 (+ getprosol)**

Pulse programme: **awhsqcedetdp-135**

Type **eda** (enter) and review the following default parameters

SW ^1H = 12 ppm, **SW ^{15}N** = **30** ppm (or other suitable values).

TD ^1H = 2K, **TD ^{15}N** = 128-160 (your choice).

NS = multiple of 8 or 16, **DS** = 16.

Type **ased** (enter) and review other parameters used in the job.

O1P = ^1H spectral window midpoint = 4.7 ppm other value of your choice.

O2P = ^{15}N spectral window midpoint = 118 ppm other value of your choice.

D1 = repetition delay = **1.0 sec** or other time of your choice.

CNST4 = $^1\text{J } ^{15}\text{N}$ - ^1H coupling constant = **90 Hz** or other value of your choice.

D24 = **D26** = $1/4\text{J}$ are auto calculated from **CNST4**

ZGOPTNS = Not used.

Check gradient settings.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K SI(F1) = 512 or 1K**

WDW(F1) = WDW(F2) = QSINE

SSB(F2) = SSB(F1) = 2

xfb, abs1 and abs2



File 104 ^1H - ^{15}N f2HSQC spectrum of a peptide which that has 9 amino acid units

2.12 ^1H - ^{15}N BESTF3HSQC spectrum

Parameter sets: **awbestf3hsqc** (+ **getprosol**)

Pulse programme: **b_hsqcetf3gpsi**

Type **eda** (enter) and review the following default parameters

SW ^1H = 12 ppm, **SW** ^{15}N = **30** ppm (or other suitable values).

TD ^1H = 2K, **TD** ^{15}N = 128-160 (your choice).

NS = multiple of 8 or 16, **DS** = 16.

Type **ased** (enter) and review other parameters used in the job.

O1P = ^1H spectral window midpoint = 4.7 ppm other value of your choice.

O3P = ^{15}N spectral window midpoint = 118 ppm other value of your choice.

D1 = repetition delay = **0.xx** sec or other time of your choice.

CNST4 = ^1J ^{15}N - ^1H coupling constant = **90 Hz** or other value of your choice.

ADD best set up info

ZGOPTNS = Not used.

Check gradient settings.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K SI(F1) = 512 or 1K**

WDW(F1) = WDW(F2) = QSINE

SSB(F2) = SSB(F1) = 2

xfb, abs1 and abs2



File 18 Best f3hsqc expt

2.13 ^1H - ^{15}N BESTf3TROSY spectrum

Parameter sets: **awbestf3trosy** (+ **getprosol**)

Pulse programme: **b_rosytf3gpsi**

Type **eda** (enter) and review the following default parameters

SW ^1H = 12 ppm, **SW** ^{15}N = **30** ppm (or other suitable values).

TD ^1H = 2K, **TD** ^{15}N = 128-160 (your choice).

NS = multiple of 8 or 16, **DS** = 16.

Type **ased** (enter) and review other parameters used in the job.

O1P = ^1H spectral window midpoint = 4.7 ppm other value of your choice.

O3P = ^{15}N spectral window midpoint = 118 ppm other value of your choice.

D1 = repetition delay = **0.xx** sec or other time of your choice.

CNST4 = 1J ^{15}N - ^1H coupling constant = **90 Hz** or other value of your choice.

ADD besttrosy set up info

ZGOPTNS = Not used.

Check gradient settings.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K SI(F1) = 512 or 1K**

WDW(F1) = WDW(F2) = QSINE

SSB(F2) = SSB(F1) = 2

xfb, abs1 and abs2




File 6 Best f3trosy expt

Appendixor standalone f2/f3 1H-15N set up info file


Spectrometer routings

The required spectrometer routing for a ^{15}N **f3** channel experiment is shown in Figure 1. The routing can be checked by typing: **edsp (return)**. If a red error line appears between the second (lower) **X 500 W Amplifier** and the ^{15}N **Preamplifier** can be removed by selecting the **Default** and **Save and Close** buttons, after which **O1**, **O2** if used and **O3** must be re-entered since they will have been reset to **0 Hz**. The correct routing for an **f2** channel ^{15}N experiment is shown on the next page.



Channel routings for an **f3** ^{15}N experiment.

The **f2** channel can optionally be set to **OFF** if ^{13}C pulses are not used.



Channel routings for a ^{12}C ^{15}N experiment.

Processing options

A large residual **H₂O/HOD** line when spectra are run in **9:1 H₂O:D₂O** or other solvents can be suppressed by setting up the experiment with its **01** value in **Hz**, or **01P ppm** value (typically ~ 4.7 ppm), to that of the **HOD** line and processing the spectrum using the following **ProcPars (edp)** settings:

BC_MOD = qfil

BCFW = 0.5 ppm or other suppression band width value of your choice

STSI = 0 (not used) or **1024** points when a data set acquired and processed with

2048 ¹H points to ONLY display spectral data to the left hand side (higher ppm side) of the noisy residual H₂O/HOD line.

Baseline correction			
ABSG	5	5	Degree of polynomial for abs (0..5)
ABSF1 [ppm]	1000.00000	1000.00000	Left limit for absf
ABSF2 [ppm]	-1000.00000	-1000.00000	Right limit for absf, abs1, abs2
BCFW [ppm]	0.50000	1.00000	Filter width for bc (sfil/qfil)
COROFFS [Hz]	0	0	Correction offset for BC_MOD=spol etc.
BC_mod	qfil	no	Fid baseline modes for em, ft, xfb,...

Qfil settings.

The default values are: **BC_MOD = no**, **BCFW = 0 or 1.00000**

Provided the s/n ratio of Fourier transformed ¹⁵N axis data points is reasonable, **linear prediction** can be applied to improve the apparent resolution of that axis.

Fourier transform			
TDeff	0	0	Number of fid data points used by ft
STSR	0	0	First output point of strip transform
STSI	1024	0	Total number of output points of strip transform
ME_mod	no	LPfc	Linear prediction for ft, xfb, ...
NCOEF	0	32	Number of LP coefficients
LPBIN	0	256	Number of output points for LP
TDoff	0	0	Number of back-predicted points

STSI and linear prediction settings

The default values for these parameters when they are not used are:

ME_MOD = no, **NCOEF = 0**, **LPBIN = 0**