

KJM 9250



AVIIIHD-800 MHz SELHSQC, SELHSQC-DIPSII2
and SELHMBC Experiments

Version 1.0

Topspin 3.5

Windows 7



© Professor Emeritus Alistair Lawrence Wilkins,
University of Waikato, New Zealand.
January 2018

SELHSQC + SELHMBC Experiments on the AVIIIHD-800

1.1 Introduction

aw coded 1D-SELHSQC, SELHSQC-DIPSII2 and SELHMBC experiments and parameter sets have been set up on the AVIIIHD-800 spectrometer.

1.2 NS x TD0 option

SELHSQC and SELHMBC experiments can be run using the NS x TD0 option where NS is a multiple of 8, 16, 32, 64 (etc) and TD0 is any number >1.

Do NOT use the TR command as a NS x TD0 experiment proceeds. Multiples of NS scans will be automatically saved and can be processed using the FT or EFP commands the experiment proceeds. A run can be terminated at any time using the STOP (*NB not the HALT*) command. This will ensure a multiple of 4 or 8 scans is saved as is required by some of the selective excitation experiments.

1.3 Signals with short T₂'s

¹³C signals with short T₂'s may give no result in SELHMBCCCT or SELHSQC experiments run with a 180 degree p36 = 40000 usec Q3 pulse in which case (*after* using the getprosol command) one can try *halving* the p36 pulse time to 20000 usec *and subtracting 6 db* from its sp26(db) power level.

SELHMBCCQ5 experiments are run with a parameter set saved power level for its 20000 usec 90 degree Q5 pulses which is *not* over written or updated by the getprosol command.

1.4 SELHSQC and SELHMBC Experiments

The following 1D-Selective experiments have been set up on the AVIIIHD-800 spectrometer.

2.1 SELHSQC and SELHSQCND spectra

2.2 SELHSQC-DIPSII2 and SELHSQCND-DIPSII2 spectra

2.3 SELHMBCCCT spectra

2.4 SELHMBCCQ5 spectra

2.1 SELHSQC and SELHSQCND Spectra

Parameter set: **awselhsqc** or **awselhsqcnd** (+ **getprosol**)

Pulse programme: **awselhsqcgpsisp** or **awselhsqcndgpsisp**

Prior to running a **SELHSQC** experiment run a standard ^{13}C or **DEPT** experiment and determine the **O1** frequency of the ^{13}C signal in **Hz** to be selectively excited. Enter this value as **O2 (Hz)**.

TD = **SI** = 64 K.

SW = 14 ppm, **O1P** = 6 ppm. Adjust **SW** and **O1P** as required.

O2 = frequency of the ^{13}C signal in **Hz** to be selectively excited.

NS = multiple of 8 or 16, **DS** = 4 or 8.

or **NS** x **TD0** scans where **TD0** = any positive number.


D1 = 1 sec or other value of your choice.

D24 is automatically calculated from **CNST2** ($^1J_{\text{C-H}}$).

CNST2 = $^1J_{\text{C-H}}$; typically 125 to 160 Hz for $\text{sp}^3\text{-sp}^2$ carbons. Furan or pyrrole ring carbons adjacent to hetero atoms will have $^1J = 200\text{-}220$ Hz.

Shaped pulse types and powers are read in by the **getprosol** command.

Process with **EFP** (applies **LB**, typically use 0.3- 0.5 Hz).



7-9 ppm region of 800 MHz **SELHSQC** (*lower*) and **SELHSQCND** (*upper*) spectra determined for quinine in $\text{D}_6\text{-DMSO}$ with selective excitation of the ^{13}C signal at 131.6 ppm (**O2** = 26475 Hz). The 1J correlated proton signal occurs at 7.95 ppm.

2.2 SELHSQC-DIPS12 and SELHSQCND-DIPS12 Spectra

Parameter set: **awselhsqc-dipsi2** or **awselhsqcnd-dipsi2** (+ **getprosol**)

Pulse programme: **awselhsqcgpdigpsisp** or **awselhsqcgpdigpndsisp**

Prior to running a **SELHSQC-DIPS12** experiment run a standard ^{13}C or **DEPT** experiment and determine the **O1** frequency of the ^{13}C signal in **Hz** to be selectively excited. Enter this value as **O2 (Hz)**.

TD = **SI** = 64 K.

SW = 14 ppm, **O1P** = 6 ppm. Adjust **SW** and **O1P** as required.

O2 = frequency of the ^{13}C signal in **Hz** to be selectively excited.

NS = multiple of 8 or 16, **DS** = 4 or 8.

or **NS** x **TD0** scans where **TD0** = any positive number.

D1 = 1 sec or other value of your choice.


D9 = 80 msec or other time of your choice (6-160 msec).

D24 is automatically calculated from **CNST2** ($^1J_{\text{C-H}}$).

CNST2 = $^1J_{\text{C-H}}$; typically 125 to 160 Hz for $\text{sp}^3\text{-sp}^2$ carbons. Furan or pyrrole ring carbons adjacent to hetero atoms will have $^1J = 200\text{-}220$ Hz.

Shaped pulse types and powers read in by the **getprosol** command

Process with **EFP** (applies **LB**, typically use 0.3 - 0.5 Hz).



7-9 ppm regions of the **SELHSQC-DIPS12** (*lower*) and **SELHSQCND-DIPS12** (*upper*) spectra determined for quinine in $\text{D}_6\text{-DMSO}$ with selective excitation of the ^{13}C signal at 131.6 ppm (**O2** = 26475 Hz). The 1J correlated proton signal occurs at 7.95 ppm.

Correlated ^1H NMR signals observed in coupled SELHSQCND-DIPS12 spectra show 1J , 2J , or nJ $^{13}\text{C}\text{-}^1\text{H}$ couplings depending on the number of bonds between the selectively excited ^{13}C signal and correlated proton signals.

2.4 SELHMBCQ5 Spectra

Parameter sets: **awselhmabcq5 (+ getprosol)**

Pulse programmes: **awselhmabcq5**

Prior to running a **SELHMMBC** experiment run a standard ^{13}C experiment and determine the **O1** frequency of the ^{13}C signal in **Hz** to be selectively excited. Enter this value as **O2 (Hz)**.

TD = **SI** = 32 K.

SW = 14 ppm, **O1P** = 6 ppm. Adjust **SW** and **O1P** as required.

O2 = frequency of the ^{13}C signal in **Hz** to be selectively excited.

NS = multiple of 8 or 16, **DS** = 4 or 8.


D1 = 1.5 sec or other value of your choice.

CNST2 = $^1J_{\text{C-H}}$; typically 125 to 160 Hz for $\text{sp}^3\text{-sp}^2$ carbons. Furan or pyrrole ring carbons adjacent to hetero atoms will have $^1J = 200\text{-}220$ Hz.

CNST13 = J_{LR} = 8 Hz or other value of your choice.

Type **ased** (return) and check gradients and other parameters are OK. Check that the experiment is set up to use a **20000 usec/18.66 db p32:sp63 Q5.1000 pulse** which is not prosol Table linked.

Process with **EFP** *and* **MC** (or **PS**) (applies **LB**, typically use 0.3- 0.5 Hz).



Lower: SELHMBCSINC spectrum of ivermectin in D₆-DMSO entered at O2 = 19286 Hz (95.9 ppm). **Upper:** SELHMBCQ3 centered at O2 = 19286 Hz .

2.4 SELHMBCCCT

Parameter set: **awselhmbcct (+ getprosol)**

Pulse programme: **awselhmbcct**

Prior to running a **SELHMMBCCT** experiment run a standard ^{13}C experiment and determine the **O1** frequency of the ^{13}C signal in **Hz** to be selectively excited. Enter this value as **O2 (Hz)**.

TD = **SI** = 32 K.

SW = 14 ppm, **O1P** = 6 ppm. Adjust **SW** and **O1P** as required.

O2 = frequency of the ^{13}C signal in **Hz** to be selectively excited.

NS = multiple of 8 or 16, **DS** = 4 or 8.


D1 = 1.5 sec or other value of your choice.

CNST6 = min. $^1J_{\text{C-H}}$ = 120 Hz, **CNST7** = max. $^1J_{\text{C-H}}$ = 170 Hz

CNST13 = J_{LR} = 8 Hz or other value of your choice

Type **ased** (return) and check gradients and other parameters are OK including the use of a prosol Table linked **40000 usec p36:sp26 Q3 pulse**.

Process with **EFP** *and* **MC** (or **PS**) (applies **LB**, typically use 0.3- 0.5 Hz).



Lower: SELHMBCCCT spectrum of ivermectin in D₆-DMSO entered at O2 =

2.5 SELHMBCT.3 Spectra

Parameter set: **awselhmbsct.3 (+ getprosol)**

Pulse programme: **awselhmbsct.3**

Prior to running a **SELHMMBCCT.3** experiment run a standard ^{13}C experiment and determine the **O1** frequency of the ^{13}C signal in **Hz** to be selectively excited. Enter this value as **O2 (Hz)**.

TD = **SI** = 32 K.

SW = 14 ppm, **O1P** = 6 ppm. Adjust **SW** and **O1P** as required.

O2 = frequency of the ^{13}C signal in **Hz** to be selectively excited.

NS = multiple of 8 or 16, **DS** = 4 or 8.

D1 = 1.5 sec or other value of your choice.


CNST6 = min. $^1J_{\text{C-H}}$ = 120 Hz, **CNST7** = max. $^1J_{\text{C-H}}$ = 170 Hz

CNST13 = J_{LR} = 8 Hz or other value of your choice.

The **Q3_surbop.1** pulse used in the **awselhmbsct.3 pp** is defined as **p32:sp57** where **p32=0.5*p36** [**p32** = 20000 usec, **p36** = 40000 usec = prosol linked value] and the parameter set saved **sp57** power level is **6 db greater** (= attenuated by a factor of 2) than that used for a 40000 usec Q3 pulse.

Type **ased** (return) and check gradients and other parameters are OK.

Process with **EFP** and **MC** (or **PS**) (applies **LB**, typically use 0.3- 0.5 Hz).



SELHMBCT.3 spectrum of ivermectin in D_6 -DMSO entered at **O2** = 19286 Hz (95.9 ppm).

SELHMBCT.3 spectra may include low levels of correlations ex nearby carbon signals.