



KJM 9250

^1H NMR spectra on the AVI-600 spectrometer.

Version 7.3

Topspin 1.3 Windows XP AVI600



© Professor Emeritus Alistair Lawrence Wilkins,
University of Waikato, New Zealand.
January 2020

¹H NMR spectra on the AVI-600

1.0 Introduction

aw coded ¹H NMR parameter files generally use a 90° pulse for maximum signal.

Best ¹H resolution is obtained using **FT** and **PK** (or **APK**) processing. **FT** processing does not apply a line broadening factor. **EF** or **EFP** processing applies a line broadening factor (**LB**).

Resolution enhancement uses negative **LB** values. Try **LB** = -1.0 to -2.5 Hz with **GB** = 0.33, and **GFP** processing. Remember to reset **LB** and **GB** to their normal values (0.1 and 0 respectively) after **GFP** processing.

1.1 Presaturation Experiments

Continuous wave (**CW**), pulsed presaturation (**PS**) or excitation sculptured (**ES**) can be used to presaturate ¹H NMR signals.

CW presaturation power levels (db settings) can be increased or decreased by subtracting or adding 3-12 db respectively. 6 db = a factor of 2. Pulsed presaturation (**PS**) uses a looped **P18 squa100.1000** pulse. The AVI's **prosol** table linked **P18** pulse time (**10 msec**) is different from that used on the AVII-600 or AVIIHD-800.

AVI-600 **ES** pulses are defined as **2000 usec p12:sp1 or p40:sp10 squa100.1000** pulses depending on which **prosol relations** option is used in a pulse program.

The **ES** shaped pulse's excitation window can be decreased by doubling its shaped pulse time from 2000 usec to 4000 usec and halving its power by adding 6 db to that read in using the **getprosol** command.

Bruker sometimes uses the **NOESYPR1D** pulse programme with a short **d8** time to acquire **QNMR** spectra.

2.0 ¹H NMR experiments

2.1 ¹H NMR with a 30, 45 or 90 degree pulse

2.2 ¹H NMR with CW presaturation

2.3 ¹H NMR with F1 and F2 CW presaturation

2.4 ¹H NMR with ES peak suppression

2.5 ¹H NMR with combined ES + CW presaturation on F1

2.6 ¹H NMR with combined ES + CW presaturation on F1 and CW presaturation on F2

2.7 ¹H NMR with three peak ES + dual CW presaturation

2.8 ¹H NOESYPR1D

2.9 ¹H NMR with PS (pulsed) presaturation

2.1 ^1H NMR spectra with a 30, 45 or 90 degree pulse

Parameter sets: **awproton30**, **awproton45**, **awproton90** (+ **getprosol**)

Pulse programmes: **zg30**, **awzg45** or **zg** respectively

TD = 64 K, **SI** = 64 K.

SW = 16 ppm, **O1P** = 7.0 ppm.


D1 = 1.5 sec or other time of your choice.

NS = any number, **DS** = 2, 4 or 8.

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important!*).

Process with **FT** (no line broadening) or **EFP** (applies **LB**).



AVI-600 ^1H NMR spectrum of quinine in $\text{D}_6\text{-DMSO}$.

2.2 ^1H NMR spectrum with CW presaturation

Parameter set: **awprotonpr (+ getprosol)**

Pulse programme: **zgpr**

TD = 64K, **SI** = 64 K.

SW = 18 ppm.

O1 = frequency in Hz of the F1 signal to be presaturated.
= spectral window midpoint. Check **SW** is wide enough.

PL9 = F1 presaturation power applied during **D1**.

D1 = 2 sec or other time of your choice.


Type **ased** (enter) and review parameters used in the job.

Add (or subtract) 3-12 db to **PL9** to decrease (or increase) the presaturation power.

6 db = a factor of 2. A larger attenuation setting decreases the power level.

Set **receiver gain** using **RGA** (*important!*).

Process with **EFP** (applies **LB**).



Lower: AVI-600 ^1H NMR spectrum of quinine in D_6 -DMSO.

Upper: ^1H NMR spectrum with CW presaturation of the HOD line at 3.37 ppm.

2.3 ^1H NMR spectrum with dual CW presaturation

Parameter set: **awprotonprf1prf2 (+ getprosol)**

Pulse programme: **awprotonprf1prf2**

TD = 64 K, **SI** = 64 K.

SW = 18 ppm.

O1 = frequency in Hz of the F1 signal to be presaturated.
= spectral window midpoint. Check **SW** is wide enough.

O2 = frequency in Hz of the F2 signal to be presaturated.

D1 = 2 sec or other time of your choice.

PL9 = F1 presaturation power applied during D1.


PL21 = F2 presaturation power applied during D1.

Add (or subtract) 3-12 db to **PL9** and/or **PL21** to decrease (or increase) the presaturation power. 6 db = a factor of 2. A larger attenuation setting decreases the power level.

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important!*).

Process with **EFP** (applies **LB**).



Lower: AVI-600 ^1H NMR spectrum of quinine in D_6 -DMSO.

Upper: ^1H NMR spectrum with CW presaturation of the HOD (3.37 ppm) and DMSO (2.5 ppm) lines.

2.4 ^1H NMR spectrum with ES peak suppression

Parameter set: **awprotones (+ getprosol)**

Pulse programme: **zgesgp**

TD = 64 K, **SI** = 64 K.

SW = 18 ppm.

O1 = frequency in Hz of the F1 signal to be ES suppressed.

= spectral window midpoint. Check. **SW** is wide enough.


D1 = 1.5 sec or other time of your choice.

Type **ased** (enter) and review parameters used in the job. Verify gradients are OK.

Check **P12 = 2000 usec**, **SPNAM1 = squa100.1000**

Set **receiver gain** using **RGA** (*important!*).

Process with **EFP** (applies **LB**).



Lower: AVI-600 ^1H NMR spectrum of quinine in D_6 -DMSO.

Upper: ^1H NMR spectrum with ES suppression of the HOD line at 3.37 ppm.

2.5 ^1H NMR with combined ES and CW presaturation on F1

Parameter set: **awprotonespr (+ getprosol)**

Pulse programme: **awprotonespr**

TD = 64 K, **SI** = 64 K.

SW = 18 ppm.

O1 = frequency in Hz of the F1 signal to be ES suppressed.

= spectral window midpoint. Check **SW** is wide enough.

D1 = 2 sec or other time of your choice.


PL9 = F1 presaturation power applied during **D1**.

Type **ased** (enter) and review parameters used in the job. Verify gradients are OK.

Check **P12** = 2000 usec, **SPNAM1** = **squa100.1000**.

Set **receiver gain** using **RGA** (*important!*).

Process with **EFP** (applies **LB**).



Lower: AVI-600 ^1H NMR spectrum of quinine in D_6 -DMSO.

Upper: ^1H NMR spectrum with combined ES and CW presaturation of the HOD line at 3.37 ppm.

2.6 ^1H NMR spectrum with combined ES+CW presaturation on F1 and CW presaturation on F2

Parameter set: **awprotonesprf1prf2 (+ getprosol)**

Pulse programme: **awprotonesprf1prf2**

TD = 64 K, **SI** = 64 K.

SW = 18 ppm.

O1 = frequency in Hz of the F1 signal to be combined ES + CW suppressed.

= spectral window mid-point. Check SW is wide enough.

O2 = frequency in Hz of the F2 signal to be CW presaturated.

D1 = 2 sec or other time of your choice.

PL9 = F1 presaturation power applied during **D1**.


PL21 = F2 presaturation power applied during **D1**.

Type **ased** (enter) and review parameters used in the job. Verify gradients are OK.

Check **P40** = 2000 usec, **SPNAM10** = **squa100.1000**.

Set **receiver gain** using **RGA** (*important!*).

Process with **EFP** (applies **LB**).



Lower: AVI-600 ^1H NMR spectrum of quinine in D_6 -DMSO.

Upper: ^1H NMR with combined ES + CW presaturation of the HOD line (3.37 ppm) on F1 and the DMSO line (2.5 ppm) on F2.

2.7 ^1H NMR spectrum with three peak ES+ dual CW presaturation

Parameter set: **awprotonesprf1prf2 (+ getprosol)**

Pulse programme **awprotonesprf1prf2**

TD = 64 K, **SI** = 64 K.

SW = 20 ppm.

O1 = frequency in Hz of the F1 signal to be CW suppressed.

= spectral window midpoint. Check SW is wide enough.

O1* = frequency in Hz of the F1 signal to be ES suppressed

SPOFFS10 = (**O1***-**O1**) Hz (may be a positive or negative value).

O2 = frequency in Hz of the F2 signal to be CW presaturated.

D1 = 2 sec or other time of your choice.

PL9 = F1 presaturation power applied during **D1**.


PL21 = F2 presaturation power applied during **D1**.

Type **ased** (enter) and review parameters used in the job. Verify gradients are OK.

Check **P40** = 2000 usec, **SPNAM10** = **squa100.1000**.

Set **receiver gain** using **RGA** (*important!*).

Process with **EFP** (applies **LB**).



Lower: AVI-600 ^1H NMR spectrum of quinine in D_6 -DMSO.

Upper: ^1H NMR with CW presaturation on F1 of quinine's OCH_3 signal (3.89 ppm), offset ES suppression of the HOD line (3.37 ppm) and CW presaturation on F2 of the DMSO signal (2.5 ppm).

2.8 NOESYPR1D with CW presaturation

Parameter set: **awnoesypr1d (+ getprosol)**

Pulse programme: **awnoesypr1d**

TD = 32 K or 64K , **SI** = 32 K or 64 K.

SW = 18 ppm.

O1 = frequency in Hz of the F1 signal to be presaturated.

= spectral window midpoint. Check **SW** is wide enough.

D1 = 2 sec or other time of your choice.

D8 = 0.05 sec (NOESY delay) or other time of your choice.

PL9 = **F1** presaturation power applied during **D1**.


Add (or subtract) 3-12 db to **PL9** to decrease (or increase) the presaturation power.

6 db = a factor of 2. A larger attenuation setting decreases the power level.

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important!*).

Process with **EFP** (applies **LB**).



Lower: AVI-600 ¹H NMR spectrum of quinine in D₆-DMSO.

Upper: ¹H NOESYPR1D spectrum with CW presaturation of the HOD line at 3.37 ppm.

2.9 ^1H NMR spectrum with PS presaturation

Parameter set: **awprotonps (+ getprosol)**

Pulse programme: **awprotonps**

TD = 64K, **SI** = 64 K.

SW = 18 ppm.

O1 = frequency in Hz of the F1 signal to be presaturated.

= spectral window midpoint. Check **SW** is wide enough.

P18 = 10 ms (=10000 us) pulsed presaturation power applied during **D1**.

A different **P18** time is used on the AVI-600 and the AVIIIHD-800.

D1 = 2 sec or other time of your choice.


L6 = number of pulsed presaturation loops is automatically calculated from **P18** and **D1**.

SPNAM6 = presaturation pulse type = **Squa100.1000**.

Type **ased** (enter) and review parameters used in the job.

Set **receiver gain** using **RGA** (*important!*).

Process with **EFP** (applies **LB**).



Lower: AVI-600 ^1H NMR spectrum of quinine in D_6 -DMSO.

Upper: ^1H NMR spectrum with pulsed presaturation of the HOD line at 3.37 ppm.