



KJM 9250

**SHSQC, SHMBC spectra on the AVI600 Spectrometer**

Version 7.3

Topspin 1.3 Windows XP AVI600



© Professor Emeritus Alistair Lawrence Wilkins,  
University of Waikato, New Zealand.  
January 2020

# AVI-600 SHSQC and SHMBC Experiments

## 1.0 Introduction

aw coded AVI-600 SHSQC and SHMBC parameter sets are set up with 2K acquired  $^1\text{H}$  points and 128  $^{13}\text{C}$  increments.  $^1\text{H}$  and  $^{13}\text{C}$  spectral windows and their mid points must be determined in a standard  $^{13}\text{C}$  or HSQC/HMBC experiment before setting up an SHSQC or SHMBC experiments.

NB: Some of the shaped pulse powers that are used in AVI-600/TS2.1 SHSQC and SHMBC experiments *are not prosol Table linked so they are saved in parameter sets.*

## Processing

SHSQC45 and SHSQC135 experiments are phase sensitive experiments. These spectra should be phased **before** using the **abs1** and **abs2** commands. Low level  $^2J$  correlations may be observed in SHSQC spectra. SHMBCQ5 spectra are processed with **xfb** *and* **xf2m**.

The  $^{13}\text{C}$  axis resolution of SHSQC and SHMBC spectra acquired using 128 increments across a 30 ppm window is typically 4-5 times greater than that of standard full window HSQC and HMBC spectra acquired with 160-256 or more increments.

## 2.0 SHSQC, SHMBC Experiments

### 2.1 SHMBC45 spectrum

### 2.2 SHSQC135 spectrum

### 2.3 SHMBCQ5 spectrum

## 2.1 SHSQC45 Spectrum

Parameter set: **awshsqc45 (+ getprosol)**

Pulse programme: **awshsqcetgpsisp2.2-45**

**d24** is automatically calculated from **cnst2**

**SW(<sup>13</sup>C) = 30 ppm**, excitation band width is ca **20 ppm**

**TD(F2) = 2K**, **TD(F1) = 128** points or other value of your choice (64-256 points).

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **O1P = <sup>1</sup>H** spectral window midpoint in ppm.

Enter **O2P = <sup>13</sup>C** spectral window midpoint in ppm.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **1.5 sec** or other time of your choice.

**CNST2** = <sup>1</sup>J coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a **2000 usec p33:sp28 Q3.1000** pulse with a parameter set saved **PL28** power level of **~ 20 db**.


Set receiver gain using **RGA** (*Important!*).

Process with: **SI(F2) = 2K**, **SI(F1) = 512 or 1024 points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb, abs1 and abs2**



AVI-600 SHSQC spectrum of quinine in D<sub>6</sub>-DSMO. The <sup>13</sup>C axis is centered at 40 ppm.

## 2.2 SHSQC135

Parameter set: **awshsqc135 (+ getprosol)**

Pulse programme: **awshsqcedetgpsisp2.3-135**

**d21** and **d24** are automatically calculated from **cnst2**

**SW(<sup>13</sup>C) = 30 ppm**, excitation band width is ca **20 ppm**

**TD(F2) = 2K**, **TD(F1) = 128** points or other value of your choice (64-256 points).

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **O1P = <sup>1</sup>H** spectral window midpoint in pp

Enter **O2P = <sup>13</sup>C** spectral window midpoint in ppm.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **1.5 sec** or other time of your choice.

**CNST2** = <sup>1</sup>J coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a **2000 usec p33:sp28 Q3.1000** pulse with a parameter set saved **PL28** power level of **~ 20 db**.


Set receiver gain using **RGA** (*Important!*).

Process with: **SI(F2) = 2K**, **SI(F1) = 512** or **1024** points

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb, abs1 and abs2**



AVI-600 SHSQC135 spectrum of quinine in D<sub>6</sub>-DSMO. The <sup>13</sup>C axis is centered at 40 ppm.

## 2.3 SHMBC

Parameter set: **awshmbcq5 (+ getprosol)**

Pulse programme: **awshmbcq5**

**SW(<sup>13</sup>C) = 30 ppm**, excitation band width is ca **15-18 ppm**.

**TD(F2) = 2K**, **TD(F1) = 128** points or other value of your choice (64-256 points).

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **O1P = <sup>1</sup>H** spectral window midpoint in ppm.

Enter **O2P = <sup>13</sup>C** spectral window midpoint in ppm.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **1.5 sec** or other time of your choice.

**CNST2** = <sup>1</sup>J coupling constant = **145 Hz** or other value of your choice (eg 125-220 Hz).

**CNST13** = <sup>n</sup>J selection filter = **8 Hz** or other value of your choice (eg: 6-14 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a **2500 usec p17:sp30 Q5.1000** pulse with a parameter set saved **PL30** power level of **~ 19 db**.


Set receiver gain using **RGA** (*Important!*).

Process with: **SI(F2) = 2K**, **SI(F1) = 512 or 1024 points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb, abs1 and abs2**



AVI-600 SHMBCQ5 spectrum of quinine in D<sub>6</sub>-DMSO. The <sup>13</sup>C axis was centered at 60 ppm.