



**KJM 5250 and KJM 9250**  
**SHSQC and SHMBC Experiments**  
Version 4.1  
Topspin 4.3



© Professor Emeritus Alistair Lawrence Wilkins,  
University of Waikato, New Zealand.  
March 2024



© Professor Frode Rise, University of Oslo, Norway.  
March 2024

## **AVneo400 SHSQC and SHMBC Experiments.**

**(This manual is about narrow, medium and wider F1 axis SHSQC and SHMBC experiments which can be performed without using the “selinv” command and its complexities).**

### **1.0 Introduction**

aw coded **Neo-400 SHSQC** and **SHMBC** parameter sets are set up with **1K or 2K** acquired  $^1\text{H}$  points in F2 and 128  $^{13}\text{C}$  increments in F1.

$^1\text{H}$  and  $^{13}\text{C}$  spectral windows and their mid points should be determined before setting up **SHSQC** or **SHMBC** experiments. The **O1** frequency at which **PR** or **ES** is applied at should be determined accurately in **Hz**, rather than approximately in **ppm**.

Topspin's **getprosol** and **pulsecal** commands should be used to read in **prosol Table** stored **pulse times** and **powers** and adjust them based on your samples solvent and matrix /buffer effects. Optionally, for concentrated samples, but not low level samples, **pulsecal 13c** can be used to adjust  $^{13}\text{C}$  pulse times and powers.

**Neo-400 aw** coded **shsqc135** pp's have been derived from Topspin's **hsqcedetgpsisp2.3** pp with the addition of auto-calculation of **d21** and **d24** from **cnst2** (= the  $^1J^{13}\text{C}-^1\text{H}$  coupling constant: default value = 145 Hz) and the removal of an optional **p28** trim pulse.

*Constant time* (CT) **aw** coded **shmbcct** parameter sets are set up with **min/max  $^1J$**  coupling constants of **120 Hz** and **170 Hz** respectively and an **8 Hz  $^nJ$**  selection filter .

### **1.1 Processing**

**SHSQC45** and **SHQC135** experiments are phase sensitive experiments which should be manually phased before optionally using the **abs1** and **abs2** commands. Low level  $^2J$  correlations and some correlations outside the set up  $^{13}\text{C}$  SW (F1) range may be observed in **SHSQC** spectra.

Constant time **SHMBCCT** spectra acquired in echo-antiecho mode are transformed with **xfb** and **xf2m**. Phasing of these spectra is not required.

The  $^{13}\text{C}$  axis resolution of **SHSQC** and **SHMBC** spectra acquired using 128 increments and linearly processed with **512 or 1024**  $^{13}\text{C}$  axis points is typically 3-4 times greater than that of standard full window **HSQC** and **HMBC** spectra acquired with 160-256 or more increments.

### **2.0 SHSQC Experiments and Parameter Sets**

The following  $^1\text{H}$  detected **SHSQC** experiments and parameter sets have been set up on the **Neo-400** spectrometer.

- |            |                   |                          |
|------------|-------------------|--------------------------|
| <b>2.1</b> | <b>shsqc45</b>    | narrow window experiment |
| <b>2.2</b> | <b>shsqc135</b>   | narrow window experiment |
| <b>2.3</b> | <b>shsqc135pr</b> | with PR presaturation    |
| <b>2.4</b> | <b>shsqc135es</b> | with ES peak suppression |

|             |                       |                               |
|-------------|-----------------------|-------------------------------|
| <b>2.5</b>  | <b>shsqc135.m</b>     | medium window experiment      |
| <b>2.6</b>  | <b>shsqc135pr.m</b>   | with PR presaturation         |
| <b>2.7</b>  | <b>shsqc135es.m</b>   | with ES peak suppression      |
| <b>2.8</b>  | <b>shsqc135espr.m</b> | with ES + PR peak suppression |
| <b>2.9</b>  | <b>shsqc135.w</b>     | wider window experiment       |
| <b>2.10</b> | <b>shsqc135pr.w</b>   | with PR presaturation         |
| <b>2.11</b> | <b>shsqc135es.w</b>   | with ES peak suppression      |

### **3.0 SHMBC Experiments and Parameter Sets**

The following <sup>1</sup>H detected **SHMBC** experiments and parameter sets have been set up on the **Neo-400** spectrometer.

|             |                      |                               |
|-------------|----------------------|-------------------------------|
| <b>3.1</b>  | <b>shmbcct</b>       | narrow window experiment      |
| <b>3.2</b>  | <b>shmbcctpr</b>     | with PR presaturation         |
| <b>3.3</b>  | <b>shmbcctes</b>     | with ES peak suppression      |
| <b>3.4</b>  | <b>shmbcct.m</b>     | medium window experiment      |
| <b>3.5</b>  | <b>shmbcctpr.m</b>   | with PR presaturation         |
| <b>3.6</b>  | <b>shmbcctes.m</b>   | with ES peak suppression      |
| <b>3.7</b>  | <b>shmbcctespr.m</b> | with ES + PR peak suppression |
| <b>3.8</b>  | <b>shmbcct.w</b>     | wider window experiment       |
| <b>3.9</b>  | <b>shmbcctpr.w</b>   | with PR presaturation         |
| <b>3.10</b> | <b>shmbcctes.w</b>   | with ES peak suppression      |

### **4.0 Appendix How to offset ES from O1 in ESPR experiments**

## 2.1 SHSQC45 Spectrum

Parameter set: **awshsqc45 (+ getprosol + pulsecal)**

Pulse program: **awshsqcsetgpsisp2.2-45**

**d24** is automatically calculated from **cnst2**

**SW(<sup>13</sup>C) = 15 ppm**, excitation band width is **8-10 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1P** = <sup>1</sup>H spectral window midpoint in ppm.

Enter **O2P** = <sup>13</sup>C spectral window midpoint in ppm.

**TD(F2) = 1K or 2K**, **TD(F1) = 128 points**.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **1.5 sec** or other time of your choice.

**CNST2** = <sup>1</sup>J coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Type **ased** (enter) and review parameters used in the job.

Check gradients and shaped pulses are OK, including a prosol Table defined

**p43:sp32 Q3\_surbop.1** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K**, **SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb, manual phasing** and **abs1 + abs2**



Neo-400 SHSQC45 spectrum of quinine in D<sub>6</sub>-DSMO centered at 40 ppm.

## 2.2 SHSQC135 Spectrum

Parameter set: **awshsqc135 (+ getprosol + pulsecal)**

Pulse program: **awshsqc135**

**d21** and **d24** are automatically calculated from **cnst2**

**SW(<sup>13</sup>C) = 15 ppm**; excitation band width is **8-10 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1P** = <sup>1</sup>H spectral window midpoint in ppm.

Enter **O2P** = <sup>13</sup>C spectral window midpoint in ppm.

**TD(F2) = 1K or 2K, TD(F1) = 128 points.**

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **1.5 sec** or other time of your choice.

**CNST2** = <sup>1</sup>J coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined

**p43:sp32 Q3\_surbop.1** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K, SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb, manual phasing** and **abs1 + abs2**



Neo-400 SHSQC135 spectrum of quinine in D<sub>6</sub>-DSMO centered at 40 ppm. Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.

## 2.3 SHSQC135pr Spectrum

Parameter set: **awshsqc135pr (+ getprosol + pulsecal)**

Pulse program: **awshsqc135pr**

**d21** and **d24** are automatically calculated from **cnst2**

**SW(<sup>13</sup>C) = 15 ppm**, excitation band width is **8-10 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1 = <sup>1</sup>H** spectral window midpoint in **Hz** (for **PR**).

Enter **O2P = <sup>13</sup>C** spectral window midpoint in ppm.

**TD(F2) = 1K or 2K**, **TD(F1) = 128 points**.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **2 sec** or other time of your choice.

**PLW9(db) = PR** power applied during **D1**. If required the **PR** power can be increased by subtracting 6 or 12 db from its prosol Table value.

**CNST2 = <sup>1</sup>J** coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Check that gradients and shaped pulses are OK, including a prosol Table defined

**p43:sp32 Q3\_surbop.1** pulse


Set receiver gain using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K**, **SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb, manual phasing and abs1 + abs2**



**Neo-400 SHSQC135pr** spectrum of quinine in D<sub>6</sub>-DSMO centered at 40 ppm.

Presaturation was applied at the HOD line frequency (3.35 ppm). Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.

## 2.4 SHSQC135es Spectrum

Parameter set: **awshsqc135es (+ getprosol + pulsecal)**

Pulse program: **awshsqc135es**

**d21** and **d24** are automatically calculated from **cnst2**

**SW(<sup>13</sup>C) = 15 ppm**, excitation band width is **8-10 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1 = <sup>1</sup>H** spectral window midpoint in **Hz** (for **ES**).

Enter **O2P = <sup>13</sup>C** spectral window midpoint in ppm.

**TD(F2) = 1K or 2K**, **TD(F1) = 128 points**.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **1.5 sec** or other time of your choice.

**CNST2** = <sup>1</sup>J coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined

**2000 usec p40:sp10 Sinc1.1000 ES** pulse and a **p43:sp32 Q3\_surbop.1** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K**, **SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb, manual phasing** and **abs1 + abs2**



Neo-400 SHSQC135es spectrum of quinine in D<sub>6</sub>-DSMO centered at 40 ppm. ES was applied at 3.35 ppm (= the HOD line). <sup>1</sup>H signals /correlations located 0.5-0.7 ppm either side of the ES suppressed line have reduced intensity. Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.

## 2.5 SHSQC135.m Spectrum

Parameter set: **awshsqc135.m (+ getprosol + pulsecal)**

Pulse program: **awshsqc135.m**

**d21** and **d24** are automatically calculated from **cnst2**

**SW(<sup>13</sup>C) = 40 ppm**; excitation band width is ~ **30-34 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1P** = <sup>1</sup>H spectral window midpoint in ppm.

Enter **O2P** = <sup>13</sup>C spectral window midpoint in ppm.

**TD(F2) = 1K or 2K, TD(F1) = 128 points.**

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **1.5 sec** or other time of your choice.

**CNST2** = <sup>1</sup>J coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined **p24:sp16 Q3\_surbop.1** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K, SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb, manual phasing** and **abs1 + abs2**



**Neo-400 SHSQC135.m** spectrum of quinine in D<sub>6</sub>-DSMO centered at 40 ppm.

Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.



## 2.6 SHSQC135pr.m Spectrum

Parameter set: **awshsqc135pr.m (+ getprosol + pulsecal)**

Pulse program: **awshsqc135pr.m**

**d21** and **d24** are automatically calculated from **cnst2**

**SW(<sup>13</sup>C) = 40 ppm**; excitation band width is **30-34 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1 = <sup>1</sup>H** spectral window midpoint in **Hz** (for **PR**).

Enter **O2P = <sup>13</sup>C** spectral window midpoint in ppm.

**TD(F2) = 1K or 2K**, **TD(F1) = 128 points**.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **2 sec** or other time of your choice.

**PLW9(db) = PR** power applied during **D1**. If required the **PR** power can be increased by subtracting 6 or 12 db from its prosol Table value.

**CNST2 = <sup>1</sup>J** coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined

**p24:sp16 Q3\_surbop.1** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K**, **SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb, manual phasing** and **abs1 + abs2**



**Neo-400 SHSQC135pr.m** spectrum of quinine in D<sub>6</sub>-DSMO centered at 40 ppm.

Presaturation was applied at the HOD line frequency (3.35 ppm). Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.

## 2.7 SHSQC135es.m Spectrum

Parameter set: **awshsqc135es.m (+ getprosol + pulsecal)**

Pulse program: **awshsqc135es.m**

**d21** and **d24** are automatically calculated from **cnst2**

**SW(<sup>13</sup>C) = 40 ppm**; excitation band width is **30-34 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1 = <sup>1</sup>H** spectral window midpoint in **Hz** (for **ES**).

Enter **O2P = <sup>13</sup>C** spectral window midpoint in ppm.

**TD(F2) = 1K or 2K**, **TD(F1) = 128 points**.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **1.5 sec** or other time of your choice.

**CNST2** = <sup>1</sup>J coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined

**2000 usec p12:sp1 Sinc1.1000 ES** pulse and a **p24:sp16 Q3\_surbop.1** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K**, **SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb, manual phasing** and **abs1 + abs2**



Neo-400 SHSQC135es.m spectrum of quinine in D<sub>6</sub>-DSMO centered at 40 ppm. **ES** was applied at 3.35 ppm (= the HOD line). <sup>1</sup>H signals /correlations located 0.5-0.7 ppm either side of the **ES** suppressed line have reduced intensity. Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.

## 2.8 SHSQ135espr.m Spectrum

Parameter set: **awshsq135espr.m** (+ **getprosol** + **pulsecal**)

Pulse program: **awshsq135espr.m**

**d21** and **d24** are automatically calculated from **cnst2**

**SW**(<sup>13</sup>C) = **40 ppm**; excitation band width is **30-34 ppm**

Type **eda** (enter) and enter **SW**(<sup>1</sup>H) and **SW** (<sup>13</sup>C) in ppm.

Enter **O1** = <sup>1</sup>H spectral window midpoint in **Hz** (for **ES** + **PR**).

**ES** can optionally be offset from **PR** as described in the Appendix.

Enter **O2P** = <sup>13</sup>C spectral window midpoint in ppm.

**TD**(F2) = **1K** or **2K**, **TD**(F1) = **128** points.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **2 sec** or other time of your choice.

**CNST2** = <sup>1</sup>J coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

**PLW9**(db) = **PR** power applied during **D1**. If required the **PR** power can be increased by subtracting 6 or 12 db from its prosol Table value.

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined

**2000 usec p12:sp1 Sinc1.1000 ES** pulse and a **p24:sp16 Q3\_surbop.1** pulse.

**ES** can optionally be offset from **PR** at **O1** as described in the Appendix.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI**(F2) = **1K** or **2K**, **SI**(F1) = **512** or **1K** points

**WDW**(F1) = **WDW**(F2) = **QSINE**

**SSB**(F2) = **SSB**(F1) = **2**

**xfb**, **manual phasing** and **abs1** + **abs2**



**Neo-400 SHSQ135espr.m** spectrum of quinine in D<sub>6</sub>-DSMO centered at 40 ppm. **ES** + **PR** was applied at 3.35 ppm (= the HOD line). <sup>1</sup>H signals /correlations located 0.5-0.7 ppm either side of the **ES** suppressed line have reduced intensity. Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.

## 2.9 SHSQC135.w Spectrum

Parameter set: **awshsqc135.w (+ getprosol + pulsecal)**

Pulse program: **awshsqc135.w**

**d21** and **d24** are automatically calculated from **cnst2**

**SW(<sup>13</sup>C) = 60 ppm**; excitation band width is **45-50 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1P** = <sup>1</sup>H spectral window midpoint in ppm.

Enter **O2P** = <sup>13</sup>C spectral window midpoint in ppm.

**TD(F2) = 1K or 2K, TD(F1) = 128 points.**

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **1.5 sec** or other time of your choice.

**CNST2** = <sup>1</sup>J coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined

**p33:sp23 Q3\_surbop.1** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K, SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb, manual phasing** and **abs1 + abs2**



Neo-400 SHSQC135.w spectrum of quinine in D<sub>6</sub>-DSMO centered at 40 ppm.

Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.

## 2.10 SHSQC135pr.w Spectrum

Parameter set: **awshsqc135pr.w (+ getprosol + pulsecal)**

Pulse program: **awshsqc135pr.w**

**d21** and **d24** are automatically calculated from **cnst2**

**SW(<sup>13</sup>C) = 60 ppm**; excitation band width is **45-50 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1 = <sup>1</sup>H** spectral window midpoint in **Hz** (for **PR**).

Enter **O2P = <sup>13</sup>C** spectral window midpoint in ppm.

**TD(F2) = 1K or 2K, TD(F1) = 128 points.**

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **2 sec** or other time of your choice.

**PLW9(db) = PR** power applied during **D1**. If required the **PR** power can be increased by subtracting 6 or 12 db from its prosol Table value.

**CNST2 = <sup>1</sup>J** coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Check that gradients and shaped pulses are OK, including a prosol Table defined **p33:sp323** Q3 type pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K, SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb, manual phasing and abs1 + abs2**



Neo-400 SHSQC135pr.w spectrum of quinine in D<sub>6</sub>-DSMO centered at 40 ppm.

Presaturation was applied at the HOD line frequency (3.35 ppm). Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.

## 2.11 SHSQC135es.w Spectrum

Parameter set: **awshsqc135esw (+ getprosol + pulsecal)**

Pulse program: **awshsqc135es.w**

**d21** and **d24** are automatically calculated from **cnst2**

**SW(<sup>13</sup>C) = 60 ppm**; excitation band width is **45-50 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1 = <sup>1</sup>H** spectral window midpoint in **Hz** (for **ES**).

Enter **O2P = <sup>13</sup>C** spectral window midpoint in ppm.

**TD(F2) = 1K or 2K, TD(F1) = 128 points.**

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **1.5 sec** or other time of your choice.

**CNST2** = <sup>1</sup>J coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined

**2000 usec p12:sp1 Sinc1.1000 ES** pulse and a **p33:sp23 Q3\_surbop.1** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K, SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb, manual phasing** and **abs1 + abs2**



**Neo-400 SHSQC135es.w** spectrum of quinine in D<sub>6</sub>-DSMO centered at 40 ppm. **ES** was applied at 3.35 ppm (= the HOD line). <sup>1</sup>H signals /correlations located 0.5-0.7 ppm either side of the **ES** suppressed line have reduced intensity. Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.

### 3.1 SHMBCCT Spectrum

Parameter set: **awshmbcct (+ getprosol + pulsecal)**

Pulse program: **shmbcctetgpl2nd**

**SW(<sup>13</sup>C) = 15 ppm**, excitation band width is **8-10 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1P** = <sup>1</sup>H spectral window midpoint in ppm.

Enter **O2P** = <sup>13</sup>C spectral window midpoint in ppm.

**TD(F2) = 1K or 2K**, **TD(F1) = 128 points**.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **1.5 sec** or other time of your choice.

**CNST6 = 120 Hz**, **CNST7 = 170 Hz** = min/max <sup>1</sup>J coupling constants.

**CNST13** = <sup>n</sup>J selection filter = **8 Hz** or other value of your choice (eg: 6-14 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined **p43:sp32 Q3\_surbop.1** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K**, **SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb** *and* **xf2m** (*and* **abs1 + abs2**)



**Neo-400 SHMBCCT** spectrum of quinine in D<sub>6</sub>-DMSO centered at 154 ppm.

Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.

### 3.2 SHMBCCTpr Spectrum

Parameter set: **awshmbcctpr (+ getprosol + pulsecal)**

Pulse program: **awshmbcctpr**

**SW(<sup>13</sup>C) = 15 ppm**, excitation band width is **8-10 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1 = <sup>1</sup>H** spectral window midpoint in **Hz** (for **PR**).

Enter **O2P = <sup>13</sup>C** spectral window midpoint in ppm.

**TD(F2) = 1K or 2K**, **TD(F1) = 128 points**.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **2 sec** or other time of your choice.

**CNST6 = 120 Hz**, **CNST7 = 170 Hz** = min/max <sup>1</sup>J coupling constants.

**CNST13 = <sup>n</sup>J selection filter = 8 Hz** or other value of your choice (eg: 6-14 Hz).

**PLW9(db) = PR** power applied during **D1**. If required the **PR** power can be increased by subtracting 6 or 12 db from its prosol Table value.

Check that gradients and shaped pulses are OK, including a prosol Table defined **p43:sp32 Q3\_surbop.1** pulse


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K**, **SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb and xf2m (and abs1 + abs2)**



**Neo-400 SHMBCCTpr** spectrum of quinine in D<sub>6</sub>-DMSO centered at 154 ppm.

Presaturation was applied at the HOD line frequency (3.35 ppm). Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.



### 3.3 SHMBCCTes Spectrum

Parameter set: **awshmbcctes (+ getprosol + pulsecal)**

Pulse program: **awshmbcctes**

**SW(<sup>13</sup>C) = 15 ppm**, excitation band width is **8-10 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1** = <sup>1</sup>H spectral window midpoint in **Hz** (for **ES**).

Enter **O2P** = <sup>13</sup>C spectral window midpoint in ppm.

**TD(F2) = 1K or 2K**, **TD(F1) = 128 points**.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **1.5 sec** or other time of your choice.

**CNST6 = 120 Hz**, **CNST7 = 170 Hz** = min/max <sup>1</sup>J coupling constants.

**CNST13** = <sup>n</sup>J selection filter = **8 Hz** or other value of your choice (eg: 6-14 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined

**2000 usec p40:sp10 Sinc1.1000 ES** pulse and a **p43:sp32 Q3\_surbop.1** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K**, **SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb and xf2m (and abs1 + abs2)**



**Neo-400 SHMBCCTes** spectrum of quinine in D<sub>6</sub>-DSMO centered at 154 .. **ES** was applied at 3.35 ppm (= the HOD line). <sup>1</sup>H signals /correlations located 0.5- 0.7 ppm either side of the **ES** suppressed line have reduced intensity. Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.

### 3.4 SHMBCCT.m Spectrum

Parameter set: **awshmbcct.m (+ getprosol + pulsecal)**

Pulse program: **awshmbcct.m**

**SW(<sup>13</sup>C) = 40 ppm**, excitation band width is **28-32 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1P** = <sup>1</sup>H spectral window midpoint in ppm.

Enter **O2P** = <sup>13</sup>C spectral window midpoint in ppm.

**TD(F2) = 1K or 2K**, **TD(F1) = 128 points**.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **1.5 sec** or other time of your choice.

**CNST6 = 120 Hz**, **CNST7 = 170 Hz** = min/max <sup>1</sup>J coupling constants.

**CNST13** = <sup>n</sup>J selection filter = **8 Hz** or other value of your choice (eg: 6-14 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined **p24:sp16 Q3\_surbop.1** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K**, **SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb and xf2m (and abs1 + abs2)**



Neo-400 **SHSQCCT.m** spectrum of quinine in D<sub>6</sub>-DSMO centered at 145 ppm.

Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.

**SHMBCCTPR** (Constant time mode SHMBC experiment)

### 3.5 SHMBCCTpr.m Spectrum

Parameter set: **awshmbcctpr.m (+ getprosol + pulsecal)**

Pulse program: **awshmbcctpr.m**

**SW(<sup>13</sup>C) = 40 ppm**, excitation band width is **30-32 ppm**.

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1 = <sup>1</sup>H** spectral window midpoint in **Hz** (for **PR**).

Enter **O2P = <sup>13</sup>C** spectral window midpoint in ppm.

**TD(F2) = 1K or 2K**, **TD(F1) = 128 points**.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **2 sec** or other time of your choice.

**CNST6 = 120 Hz**, **CNST7 = 170 Hz** = min/max <sup>1</sup>J coupling constants.

**CNST13 = <sup>n</sup>J** selection filter = **8 Hz** or other value of your choice (eg: 6-14 Hz).

**PLW9(db) = PR power** applied during **D1**. If required the **PR power** can be increased by subtracting 6 or 12 db from its prosol Table value.

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined **p24:sp16 Q3\_surbop.1** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K**, **SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb and xf2m (and abs1 + abs2)**



**Neo-400 SHMBCCTpr.m** spectrum of quinine in D<sub>6</sub>-DSMO centered at 145 ppm.

Presaturation was applied at the HOD line frequency (3.35 ppm). Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.

### 3.6 SHMBCCTes.m Spectrum

Parameter set: **awshmbcctes.m (+ getprosol + pulsecal)**

Pulse program: **awshmbcctes.m**

**SW(<sup>13</sup>C) = 40 ppm**, excitation band width is **28-32 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1 = <sup>1</sup>H** spectral window midpoint in **Hz** (for **ES**).

Enter **O2P = <sup>13</sup>C** spectral window midpoint in ppm.

**TD(F2) = 1K or 2K, TD(F1) = 128 points.**

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **1.5 sec** or other time of your choice.

**CNST6 = 120 Hz, CNST7 = 170 Hz** = min/max <sup>1</sup>J coupling constants.

**CNST13 = <sup>n</sup>J selection filter = 8 Hz** or other value of your choice (eg: 6-14 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined

**2000 usec p12:sp1 Sinc1.1000 ES pulse and a p24:sp16 Q3\_surbop.1 pulse.**


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K, SI(F1) = 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb and xf2m (and abs1 + abs2)**



Neo-400 SHMBCCTes.m spectrum of quinine in D<sub>6</sub>-DSMO centered at 145 ppm. ES was applied at 3.35 ppm (= the HOD line). <sup>1</sup>H signals /correlations located 0.5-0.7 ppm either side of the ES suppressed line have reduced intensity. Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.

### 3.7 SHMBCCTespr.m Spectrum

Parameter set: **awshmbcctwespr.w (+ getprosol + pulsecal)**

Pulse program: **awshmbcctespr.w**

**SW(<sup>13</sup>C) = 40 ppm**, excitation band width is **30-32 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1 = <sup>1</sup>H** spectral window midpoint in **Hz** (for **ES +PR**).

**ES** can optionally be offset from **PR** as described in the Appendix.

Enter **O2P = <sup>13</sup>C** spectral window midpoint in ppm.

**TD(F2) = 1K or 2K**, **TD(F1) = 128 points**.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **2 sec** or other time of your choice.

**CNST6 = 120 Hz**, **CNST7 = 170 Hz** = min/max <sup>1</sup>J coupling constants.

**CNST13 = <sup>n</sup>J** selection filter = **8 Hz** or other value of your choice (eg: 6-14 Hz).

**PLW9(db) = PR power** applied during **D1**. If required the **PR power** can be increased by subtracting 6 or 12 db from its prosol Table value.

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined

**2000 usec p12:sp1 Sinc1.1000 ES** pulse and a **p24:sp16 Q3\_surbop.1** pulse.

**ES** can optionally be offset from **PR** at **O1** as described in the Appendix.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K**, **SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb and xf2m (and abs1 + abs2)**



**Neo-400 SHMBCCTespr.m** spectrum of quinine in D<sub>6</sub>-DSMO centered at 145 ppm.

Combined **ES + PR** was applied at 3.35 ppm (= the HOD line). <sup>1</sup>H signals /correlations located 0.5-0.7 ppm either side of the **ES** suppressed line have reduced intensity.

Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity

### 3.8 SHMBCCT.w Spectrum

Parameter set: **awshmbcct.w (+ getprosol + pulsecal)**

Pulse program: **awshmbcct.w**

**SW(<sup>13</sup>C) = 60 ppm**, excitation band width is **45-50 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1P** = <sup>1</sup>H spectral window midpoint in ppm.

Enter **O2P** = <sup>13</sup>C spectral window midpoint in ppm.

**TD(F2) = 1K or 2K**, **TD(F1) = 128 points**.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **1.5 sec** or other time of your choice.

**CNST6 = 120 Hz**, **CNST7 = 170 Hz** = min/max <sup>1</sup>J coupling constants.

**CNST13** = <sup>n</sup>J selection filter = **8 Hz** or other value of your choice (eg: 6-14 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined **p33:sp23 Q3\_surbop.1** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K**, **SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb and xf2m (and abs1 + abs2)**



**Neo-400 SHMBCCT.w** spectrum of quinine in D<sub>6</sub>-DSMO centered at 140 ppm.

Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.

### 3.9 SHMBCCTpr.w Spectrum

Parameter set: **awshmbcctpr.w (+ getprosol + pulsecal)**

Pulse program: **awshmbcctpr.w**

**SW(<sup>13</sup>C) = 60 ppm**, excitation band width is **45-50 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1** = <sup>1</sup>H spectral window midpoint in **Hz** (for **PR**).

Enter **O2P** = <sup>13</sup>C spectral window midpoint in ppm.

**TD(F2) = 1K or 2K**, **TD(F1) = 128 points**.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **2 sec** or other time of your choice.

**CNST6 = 120 Hz**, **CNST7 = 170 Hz** = min/max <sup>1</sup>J coupling constants.

**CNST13** = <sup>n</sup>J selection filter = **8 Hz** or other value of your choice (eg: 6-14 Hz).

**PLW9(db) = PR power** applied during **D1**. If required the **PR power** can be increased by subtracting 6 or 12 db from its prosol Table value.

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined **p33:sp23 Q3\_surbop.1** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K**, **SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb and xf2m (and abs1 + abs2)**



**Neo-400 SHMBCCTpr.w** spectrum of quinine in D<sub>6</sub>-DSMO centered at 140 ppm.

Presaturation was applied at the HOD line frequency (3.35 ppm). Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.

### 3.10 SHMBCCTes.w Spectrum

Parameter set: **awshmbcctes.w (+ getprosol + pulsecal)**

Pulse program: **awshmbcctes.w**

**SW(<sup>13</sup>C) = 60 ppm**, excitation band width is **45-50 ppm**

Type **eda** (enter) and enter **SW(<sup>1</sup>H)** and **SW (<sup>13</sup>C)** in ppm.

Enter **O1 = <sup>1</sup>H** spectral window midpoint in **Hz** (for **ES**).

Enter **O2P = <sup>13</sup>C** spectral window midpoint in ppm.

**TD(F2) = 1K or 2K**, **TD(F1) = 128 points**.

**NS** = multiple of 4, 8 or 16, **DS** = 8 or 16.

**D1** = repetition delay = **1.5 sec** or other time of your choice.

**CNST6 = 120 Hz**, **CNST7 = 170 Hz** = min/max <sup>1</sup>J coupling constants.

**CNST13 = <sup>n</sup>J selection filter = 8 Hz** or other value of your choice (eg: 6-14 Hz).

Type **ased** (enter) and review parameters used in the job.

Check that gradients and shaped pulses are OK, including a prosol Table defined

**2000 usec p12:sp1 Sinc1.1000 ES** pulse and a **p33:sp23 Q3\_surbop.1** pulse.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 1K or 2K**, **SI(F1) = 512 or 1K points**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb and xf2m (and abs1 + abs2)**



**Neo-400 SHMBCCTes.w** spectrum of quinine in D<sub>6</sub>-DSMO centered at 140 ppm. **ES** was applied at 3.35 ppm (= the HOD line). <sup>1</sup>H signals /correlations located 0.5-0.7 ppm either side of the **ES** suppressed line have reduced intensity. Correlations at the edges of the <sup>13</sup>C spectral window have reduced intensity.



#### 4.0 How to offset ES from O1 in an ESR experiment

By default **ES** and **PR** are applied at **O1** (Hz) frequency in aw coded **HSQC** and **HMBCCT ESR** experiments. Combined (double) **ES + PR** can be used to suppress a large HOD or solvent peak.

The **ES** pulse in **shmbc135espr.m** or **shmbcctespr.m** experiments is defined as an f1 (<sup>1</sup>H) channel **2000 usec Sinc1.1000 p12:sp1** pulse.


The frequency (in Hz) at which **ES** is applied can optionally be offset from **O1** (= the frequency PR is applied) so you can suppress two solvent lines by entering an **SPOFFS1(Hz)** offset value in its **ased** visible cell.

For example if, the **HOD line** occurs at **2701 Hz** and the **DMSO line** occurs at **2007 Hz**, the offset of the **DMSO line** (*to be ES suppressed*) relative to that of the **HOD line** (*PR suppressed at O1*) is calculated as:

$$\text{SPOFFS1(Hz)} = \text{ES offset signal (Hz)} - \text{O1 frequency (Hz)}$$

$$\text{ie } 2007 \text{ Hz (DMSO)} - 2701 \text{ Hz (HOD)} = -694 \text{ Hz}$$

The offset is negative in this case since the **DMSO line** occurs at **-694 Hz less** than that of the frequency at which **PR** is applied to the **HOD line** at **O1 Hz**.



**p12:sp1** (<triple> and <triple2>) ES pulse parameters. The **SPOFFS1(Hz)** offset value is red arrowed.

*Footnote:* **ES** pulses are defined as a **p12:sp1** pulses in pp's that run with *prosol relations* = <triple> or <triple2> or as **p40:sp10** pulses with *prosol relations* = <default> or not shown/not included as a pp line.