



**KJM 5250 and KJM 9250
HSQC135 and HMBC Experiments with Excitation Sculptured
(ES) Peak Suppression on the AVneo400 spectrometer.**

Version 3.1
Topspin 4.3



© Professor Emeritus Alistair Lawrence Wilkins,
University of Waikato, New Zealand.
March 2024



© Professor Frode Rise, University of Oslo, Norway.
March 2024

AVneo400 HSQC135 and HMBC Experiments with Excitation Sculptured (ES) Peak Suppression

1.1 Spectral Window Set up

The spectral window width (in ppm), centered around the peak to be ES suppressed should be determined in a standard ^1H NMR spectrum *before* setting up a **HSQCES** or **HMBCCTES** experiment. There should be no signals within 0.5 ppm of the upper or lower limits of the spectral window. The frequency of the peak to be ES suppressed should be determined in **Hz**.

The **getprosol** command will read in **prosol Table** saved pulse time, powers and types into an experiment. After using the **getprosol** command the **pulsecal** and **pulsecal 13C** commands can be used to adjust the ^1H and ^{13}C 90 degree pulse times and **prosol Table** linked pulse powers to take account of solvent and/or buffer matrix effects which influence a samples 90 degree pulse time.

1.2 ES set up

ES is applied at **O1 Hz** = the spectral window midpoint, or it can *optionally be offset* at **O1* Hz** where **SPOFFS1** (or **SPOFFS10** in some ES pp's) = **O1*-O1 Hz**. *The use of a 2000 usec Sincl.1000 ES pulse suppresses signals $\sim \pm 0.7$ ppm (~ 280 Hz) either side of its frequency. This band width is 2-3 times greater than that of PR presaturation.*

1.3 Processing

HSQC135ES experiments are phase sensitive experiments which should be phased **before** using the **abs1** and **abs2** (and optional **syma**) commands.

The **HMBCES** experiment is acquired in MC (magnitude) mode and does not require phasing

The **HMBCCTES** (constant time) experiment is acquired in phase sensitive mode and transformed to afford an absolute value spectrum using the **xfb** and **xf2m** commands.

2.0 Experiments and Parameter Sets

2.1 HSQC135ES Spectrum

2.2 HMBCES Spectrum

2.3 HMBCCTES (constant time) Spectrum

2.7 HSQC135ES with Excitation Sculpting

Parameter set: **awhsc135es (+ getprosol)**

Pulse programme: **awhscdetgpsisp2.3es**

d21 and **d24** are automatically calculated from **cnst2**

Type **eda** (enter) and enter **SW (¹H)** and **SW (¹³C)** in ppm.

Enter **O1** in **Hz** of the signal to be **ES** suppressed

O1 = spectral window midpoint. Check **SW (¹H)** is wide enough.

Enter **O2P** = ¹³C spectral window midpoint in ppm.

TD(F2) = 1K or 2K, **TD(F1)** = 160, or 128-256 (your choice).

NS = multiple of 4, 8 or 16, **DS** = 8 or 16.

D1 = repetition delay = 1.5-2 **sec** or other time of your choice.

CNST2 = ¹J coupling constant = **145 Hz** or other value of your choice (eg: 125-160 Hz).

Type **ased** (enter) and review parameters used in the job.

Verify that a **2000 usec sinc1.1000** shaped pulse is used.

Check gradients are OK.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = SI(F1) = 1K** or **2K**

WDW(F1) = WDW(F2) = QSINE

SSB(F2) = SSB(F1) = 2

xfb, abs1 and abs2



Neo400 HSQC135ES spectrum of quinine in D₆-DMSO with the OCH₃ signal at 3.9 ppm suppressed. If the HOD signal at 3.31 ppm was **ES** suppressed correlations located ± 0.7 ppm either side of the HOD peak would also have also suppressed.

2.8 HMBCEs with Excitation Sculpting

Parameter set: **awhmbcpr (+ getprosol)**

Pulse programme: **awhmbcgplpndqfpr**

Type **eda** (enter) and enter **SW (¹H)** and **SW (¹³C)** in ppm.

Enter **O1** in **Hz** of the signal to be **ES** suppressed

O1 = spectral window midpoint. Check **SW (¹H)** is wide enough.

Enter **O2P** = ¹³C spectral window midpoint in ppm.

TD(F2) = 1K or 2K, **TD(F1)** = 160, or 128-256 (your choice).

NS = 4, 8, 16 (multiple of 4 or 8 recommended), **DS** = 8 or 16.

D1 = repetition delay = **1.5 sec** or other time of your choice.

CNST2 = ¹J coupling constant = **145 Hz** or other value of your choice.

CNST13 = ⁿJ selection filter = **8 Hz** or other value of your choice.

Type **ased** (enter) and review parameters used in the job.

Verify that a **2000 usec sinc1.1000** shaped pulse is used.

Check gradients are OK.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = SI(F1) = 1K** or **2K**

WDW(F1) = WDW(F2) = SINE

SSB(F2) = SSB(F1) = 0

xfb, abs1 and abs2



Neo400 HMBCEs spectrum of quinine in D₆-DMSO with **ES** suppression of the DMSO signal at 2.52 ppm. Correlations located ± 0.7 ppm either side of the DMSO signal are also suppressed.

2.9 HMBCCTES with Excitation Sculpting

Parameter set: **awhmbcctes (+ getprosol)**

Pulse programme: **awhmbcctes**

Type **eda** (enter) and enter **SW (¹H)** and **SW (¹³C)** in ppm.

Enter **O1 in Hz** of the signal to be ES suppressed

O1 = spectral window midpoint. Check **SW (¹H)** is wide enough.

Enter **O2P = ¹³C** spectral window midpoint in ppm.

TD(F2) = 1K or 2K, **TD(F1)** = 128-256 (your choice).

NS = 4, 8, 16 (multiple of 4 or 8 recommended), **DS** = 8 or 16.

D1 = repetition delay = **1.5 sec** or other time of your choice.

CNST6 = min ¹J coupling constant = **120 Hz** or other value of your choice.

CNST7 = max ¹J coupling constant = **170 Hz** or other value of your choice.

CNST13 = ⁿJ selection filter = **8 Hz** or other value of your choice.

Type **ased** (enter) and review parameters used in the job.

Verify that a **2000 usec sinc1.1000** shaped pulse is used.

Check gradients are OK.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = SI(F1) = 1K or 2K**

WDW(F1) = WDW(F2) = SINE

SSB(F2) = SSB(F1) = 0

xfb, xf2m, abs1 and abs2



Neo400 HMBCCTES spectrum of quinine in D₆-DMSO with **ES** suppression of the DMSO signal at 2.52 ppm. Correlations located ± 0.7 ppm either side of the DMSO signal are also suppressed.