



**KJM 5250 and KJM 9250**  
**One Dimensional SELHSQC and SELHMBC Experiments on**  
**the AVneo400 Spectrometer.  $^1\text{H}$  detected with selective  $^{13}\text{C}$**   
**excitation.**

Version 3.1  
Topspin 4.3



© Professor Emeritus Alistair Lawrence Wilkins,  
University of Waikato, New Zealand.  
March 2024



© Professor Frode Rise, University of Oslo, Norway.  
March 2024

## AVneo400 SELHSQC + SELHMBCT Experiments

### 1.1 Introduction

Prior to running SELHSQC and SELHMBCT experiments determine the frequency of the  $^{13}\text{C}$  signal to be selectively excited in Hz and enter this value as **O2 (Hz)**.

SELHSQC and SELHMBCT experiments use a selective (soft) 180 degree refocusing **Q3** type pulse.

The use of  $^{13}\text{C}$  decoupling during FID acquisition in SELHSQC experiments on the Neo-400 results in a modest degree of probe warming with some loss of shim and resolution drift unless the FID acquisition time is short (ca 0.5 sec). Accordingly the **awselhsqc** experiment is set up as a 15.6 ppm  $^1\text{H}$  window expt with 6400 data points (AQ = 0.512 sec) which when processed can be zero filled to 16-32K points (your choice) with optional linear prediction to 12-16K data points.

Decoupling related probe heating does not occur when running SELHSQCND experiments (ND = not decoupled).

The SELHMBCT experiment (which does not use  $^{13}\text{C}$  decoupling) *only affords correlations ex quaternary (s) carbons.*

### 1.2 Signals with short $T_2$ 's

$^{13}\text{C}$  signals with short  $T_2$ 's may give no result in a SELHSQC or SELHMBCT experiments run with a **180° inversion p36 = 40000 usec Q3 pulse** in which case (*after* using the getprosol command) one can try *halving* the **p36** pulse time to **20000 usec** and *subtracting 6 db* from its **sp26(db)** power level.

### 1.3 SELHSQC and SELHMBCT Experiments

The following 1D-Selective HSQC and HMBCT experiments have been set up on the Neo-400 spectrometer.

- |                      |                                           |
|----------------------|-------------------------------------------|
| <b>2.1 SELHSQC</b>   | with $^{13}\text{C}$ decoupling during AQ |
| <b>2.2 SELHSQCND</b> | without $^{13}\text{C}$ decoupling        |
| <b>2.3 SELHMBCT</b>  | quaternary carbons only                   |

## 2.1 SELHSQC Spectrum

Parameter set: **awselhsqc** (+ getprosol)

Pulse programme: **awselhsqcgpsisp**

Prior to running a **SELHSQC** experiment run a standard  $^{13}\text{C}$  or **DEPT** experiment and determine the **O1** frequency of the  $^{13}\text{C}$  signal in **Hz** to be selectively excited.

Enter this value as **O2 (Hz)** (*Important*)

**TD** = 6400, **SI** = 16 or 32 K.

Check **AQ** is  $\sim 0.5$  sec (or less) (*Important*)

**SW** = 15.6 ppm, **O1P** = 6 ppm.

**NS** = multiple of 8 or 16, **DS** = 4 or 8.


**D1** = 1.5-2 sec or other value of your choice.

**D24** is automatically calculated from **CNST2** ( $^1J_{\text{C-H}}$ ).

**CNST2** =  $^1J_{\text{C-H}}$ ; typically 125 to 160 Hz for  $\text{sp}^3\text{-sp}^2$  carbons. Furan or pyrrole ring carbons adjacent to hetero atoms will have  $^1J = 200\text{-}220$  Hz.

Shaped pulse types and powers are read in by the **getprosol** command.

Process with **EFP** (applies **LB**, typically use 0.3- 0.5 Hz).



Neo-400 **SELHSQC** spectrum of quinine in  $\text{D}_6\text{-DMSO}$  (6-10 ppm region) with selective excitation of the  $^{13}\text{C}$  signal at 131.6 ppm (**O2** = 13244.1 Hz). The  $^1J$  correlated proton signal of this carbon occurs at 7.95 ppm.

## 2.2 SELHSQCND Spectrum

Parameter set: **awselhsqnd** (+ getprosol)

Pulse programme: **awselhsqndgpsisp**

Prior to running a **SELHSQCND** experiment run a standard  $^{13}\text{C}$  or **DEPT** experiment and determine the **O1** frequency of the  $^{13}\text{C}$  signal in **Hz** to be selectively excited.

Enter this value as **O2 (Hz)** (*Important*).

**TD** = **SI** = 16 or 32 K.

**SW** = 15.6 ppm, **O1P** = 6 ppm.

**NS** = multiple of 8 or 16, **DS** = 4 or 8.


**D1** = 1.5-2 sec or other value of your choice.

**D24** is automatically calculated from **CNST2** ( $^1J_{\text{C-H}}$ ).

**CNST2** =  $^1J_{\text{C-H}}$ ; typically 125 to 160 Hz for  $\text{sp}^3\text{-sp}^2$  carbons. Furan or pyrrole ring carbons adjacent to hetero atoms will have  $^1J = 200\text{-}220$  Hz.

Shaped pulse types and powers are read in by the **getprosol** command.

Process with **EFP** (applies **LB**, typically use 0.3- 0.5 Hz).



**Lower:** Neo400 SELHSQCND spectrum of quinine in  $\text{D}_6\text{-DMSO}$  (6-10 ppm region) with selective excitation of the  $^{13}\text{C}$  signal at 131.6 ppm (**O2** = 13244.1 Hz).

**Upper:**  $^{13}\text{C}$  decoupled SELHSQC spectrum of quinine

## 2.3 SELHMCCT Spectrum

Parameter set: **awselhmbcct** (+ getprosol)

Pulse programme: **awselhmbcct**

Prior to running a **SELHMCCT** experiment run a standard  $^{13}\text{C}$  or **DEPTQ** experiment and determine the **O1** frequency of the  $^{13}\text{C}$  signal in **Hz** to be selectively excited.

Enter this value as **O2 (Hz)** (*Important*)

**TD** = **SI** = 16 or 32 K.

**SW** = 15.6 ppm, **O1P** = 6 ppm.

**NS** = multiple of 8 or 16, **DS** = 4 or 8.

**D1** = repetition delay = **1.5 sec** or other time of your choice.


**CNST6** = min  $^1J$  coupling constant = **120 Hz** or other value of your choice.

**CNST7** = max  $^1J$  coupling constant = **170 Hz** or other value of your choice.

**CNST13** =  $^nJ$  selection filter = **8 Hz** or other value of your choice.

Shaped pulse types and powers are read in by the **getprosol** command.

Process with **EFP** and **MC** (applies **LB**, typically use 0.3- 0.5 Hz).



Neo400 **SELHMCCT** spectrum of quinine in  $\text{D}_6\text{-DMSO}$  (6-10 ppm region) with selective excitation of the  $^{13}\text{C}$  signal at 157.2 ppm.