



**KJM 5250 and KJM 9250  
AVneo400  $^{15}\text{N}$  HSQC and HMBC Spectra on the AVneo400  
Spectrometer.**

Version 3.1  
Topspin 4.3



© Professor Emeritus Alistair Lawrence Wilkins,  
University of Waikato, New Zealand.  
March 2024



© Professor Frode Rise, University of Oslo, Norway.  
March 2024

## **AVneo400 <sup>15</sup>N f2 (f3 on the AVneo800) HSQC and HMBC Spectra.**

**(The f2 or f3 radio channels give rise to the F1 dimension in the 2 D spectra, the f1 radio channel give rise to the F2 dimension in the 2D spectra)**

### **1.1 Introduction**

aw coded <sup>15</sup>N f2 HSQC, HSQC-DIPSI2 and HMBC parameter sets have been set up with 2048 (2K) points acquired across a **9 or 14 ppm <sup>1</sup>H** window and a **400 ppm <sup>15</sup>N window (SW)**. These settings were used to acquire <sup>15</sup>N-HSQC and HMBC spectra ex an aromatic compound containing an NH and a quaternary N atom. Different <sup>1</sup>H and <sup>15</sup>N settings may (will) be required for other nitrogen containing compounds such as peptides.

### **1.2 Processing**

<sup>15</sup>N f2 HSQC and HSQC-DIPSI2 experiments are phase sensitive experiments which should be phased **before** using the **abs1** and **abs2** commands.

The <sup>15</sup>N f2 HMBC experiment is an absolute value experiment. Phasing is not required.

## **2.0 Experiments and Parameter Sets**

The following aw coded <sup>15</sup>N-<sup>1</sup>H f2 HSQC, HSQC-DIPSI2 and HMBC parameter sets are available on the Neo-400 MHz spectrometer.

**2.1 <sup>15</sup>N f2 HSQC**

**2.2 <sup>15</sup>N f2 HSQC-DIPSI2**

**2.3 <sup>15</sup>N f2 HMBC**

## 2.1 $^1\text{H}$ - $^{15}\text{N}$ HSQC spectrum

Parameter set: **aw15nf2hsqc (+ get prosol)**

Pulse program: **hsqcetgpsi**

**SW  $^1\text{H}$  = 14 ppm** (or other suitable value).

**SW  $^{15}\text{N}$  = 400 ppm** (or other suitable value).

**TD  $^1\text{H}$  = 2K, TD  $^{15}\text{N}$  = 128-160** (your choice).

**NS** = multiple of 8, **DS** = 8 or 16.

Type **ased** (enter) and review other parameters used in the job.

**O1P** =  $^1\text{H}$  spectral window midpoint = 6.5 ppm other value of your choice.

**O2P** =  $^{15}\text{N}$  spectral window midpoint = 150 ppm other value of your choice.

**D1** = repetition delay = 1-2 sec or other time of your choice.

**CNST2** =  $^1\text{J}^{15}\text{N}$ - $^1\text{H}$  coupling constant = **90 Hz** or other value of your choice.

**ZGOPTNS** = Not used.

Type **ased** (enter) and review parameters used in the job.

Check **D24 =  $1/8J$  (~ 1.39 msec)** for **CNST2 = 90 Hz**.

This value is not auto-calculated using Topspin's `hsqcdiedetgpsisp.2 pp`.

Check gradients and shaped pulses are OK.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K SI(F1) = 512 or 1K**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb, abs1 and abs2**



$^{15}\text{N}$  f2 HSQC spectrum of an aromatic NH containing compound.

## 2.2 $^1\text{H}$ - $^{15}\text{N}$ HSQC-DIPS12 spectrum

Parameter sets: **aw15nf2hsqc-dipsi2 (+ getprosol)**

Pulse programme: **hsqcdietgpsi**

**SW  $^1\text{H}$  = 14 ppm** (or other suitable value).

**SW  $^{15}\text{N}$  = 400 ppm** (or other suitable value).

**TD  $^1\text{H}$  = 2K, TD  $^{15}\text{N}$  = 128-160** (your choice).

**NS = multiple of 8 or 16, DS = 16.**

**O1P =  $^1\text{H}$  spectral window midpoint = 6.5 ppm** other value of your choice.

**O2P =  $^{15}\text{N}$  spectral window midpoint = 150 ppm** other value of your choice.

**D1 = repetition delay = 1.0 sec** or other time of your choice.

**CNST2 =  $^1\text{J}^{15}\text{N}$ - $^1\text{H}$  coupling constant = 90 Hz** or other value of your choice.

**ZGOPTNS = Not used.**

Type **ased** (enter) and review parameters used in the job.

Check **D24 =  $1/8\text{J}$**  ( $\sim 1.39$  msec for **CNST2 = 90 Hz**).

This values is not auto-calculated using Topspin's hsqcdietgpsi pp.

Check gradients and shaped pulses are OK.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K SI(F1) = 512 or 1K**

**WDW(F1) = WDW(F2) = QSINE**

**SSB(F2) = SSB(F1) = 2**

**xfb, abs1 and abs2**



$^{15}\text{N}$  f2 HSQC-DIPS12 spectrum of an aromatic NH containing compound.

## <sup>15</sup>N F2 HMBC spectrum

Parameter sets: **aw15nf2hmbc+ getprosol)**

Pulse programme: **hmbcfigndqf**

Type **eda** (enter) and review the following default parameters

**SW <sup>1</sup>H = 9 ppm, SW <sup>15</sup>N = 400 ppm** (or other suitable values).

**TD <sup>1</sup>H = 2K, TD <sup>15</sup>N = 128-160** (your choice).

**NS = multiple of 8 or 16, DS = 16.**

**O1P = <sup>1</sup>H spectral window midpoint = 4.79 ppm** other value of your choice.

**O2P = <sup>15</sup>N spectral window midpoint = 100 ppm** other value of your choice.

**D1 = repetition delay = 1.0 sec** or other time of your choice.

**CNST13 = 5 Hz** (long range coupling constant).

**ZGOPTNS = Not used.**

Check gradient settings.


Set **receiver gain** using **RGA** (*Important!*).

Process with: **SI(F2) = 2K SI(F1) = 512 or 1K**

**WDW(F1) = WDW(F2) = SINE**

**SSB(F2) = SSB(F1) = 0**

**xfb and abs1 + abs2**



<sup>15</sup>N f2 HMBC spectrum of an aromatic NH and quaternary N containing compound.