



KJM 9250

Proton T_1 Spectra on the AVI-600 and AVII-600

Version 1.0



© Professor Emeritus Alistair Lawrence Wilkins,
University of Waikato, New Zealand.
February 2018

Proton T₁ Spectra on the AVI-600 and AVII-600

1.0 Introduction

An **awprotont1** parameter set and linked **VDLIST** file have been set up on the AVI-600 and AVII-600 spectrometers running under TS2.1 and TS3.2 respectively.

The **awprotont1** parameter set has **D1 = 10 sec** and a linked **VDLIST** file with 8 x delays in the range 0.03 to 10 sec with the longest delay (10 sec) the first used VDLIST value and the shortest delay (0.03 sec) the last used value.

D1 and the longest (first) **VDLIST** value should be 3-5 times the longest T₁ in the sample compound.

The **awproton16t1** parameter set with **VD_{max}** and **D1 = 16 sec** can be used for signals with longer T₁'s


T₁ data sets can be processed on the spectrometer terminal or off line using any version of Topspin provided a copy of the linked **VDLIST** file is resident in the experiment's top level folder (= default set up: see below) or a copy of the originally named **VDLIST** file is recreated in the offline terminal's *C:\Bruker\Topspin...\exp\stan\nmr\lists\vd* folder.

	uxnmr.info	2/01/2017 10:50 p...	INFO File	4 KB
	uxnmr.par	2/01/2017 10:50 p...	PAR File	24 KB
	vdlist	10/01/2018 8:49 a...	File	1 KB

NB: The experimental copy of the variable delay file is named as **vdlist** irrespective of the name of the source vdlst file.

2.0 Experiment Set Up


- 1) Create an experiment and read in the **awprotont1** parameter set (+ **getprosol**).
Pulse programme = **t1ir**.
- 2) Review default settings. These settings can be adjusted if required.
TD(F2) = 32K, TD(F1) = 8.
SW = 14 ppm, O1P = 6.3 ppm.
D1 = 10 sec.
NS = multiple of 4 or 8, DS = 4 or 8.
- 3) Type **ased** (enter) and review other parameters used in the job including the linked **VDLIST** file = **AWPROTONT1-8**
- 4) The **VDLIST** file should have the entries shown below (next page). Values are in seconds. Do not alter the values in this file. An alternatively named **VDLIST** file should be created if different **VD** and **D1** values are required for a particular compound.



- 5) Set receiver gain using **RGA** (*important!*).
- 6) Type **edp** (enter) and check that **SI(F2) = 32K**, **SI(F1) = 8**, **WDW(F2) = EM**, **LB (F2) = 0.3-0.5 Hz** or other value of your choice.

	F2	F1	Frequency axis
Reference			
SI	32768	8	Size of real spectrum
SF [MHz]	600.1300000	600.1300000	Spectrometer frequency
OFFSET [ppm]	13.31307	7.29970	Low field limit of spectrum
SR [Hz]	0	0	Spectrum reference frequency
HZpPT [Hz]	0.256882	0.249924	Spectral resolution
SPECTYP	UNDEFINED		Type of spectrum e.g. COSY, H...
Window function			
WDW	EM	SINE	Window functions for trf, xfb,...
LB [Hz]	0.30	0.30	Line broadening for em

- 7) Start the acquisition using the **ZG** command.
- 8) When the experiment has run type **rser 1** (enter) to read in the first serial file which will appear in a TEMP screen display window.
- 9) Type **EFP** (return) to transform it and phase it as per a normal proton spectrum.



- 10) Type **edp** (enter) and note the phase constants for this spectrum.


Phase correction	
PHC0 [degrees]	21.768
PHC1 [degrees]	-46.451
PH_mod	pk

- 11) Close the **TEMP** window and reload the T₁ data set file.

- 12) Type **edp** (enter) and enter the phase constants noted in step 10 above into the **F2 PHC0** and **F2 PHC1** cells and check **PH_MOD = pk**. **F1** cell info is not used.


Phase correction			
PHC0 [degrees]	21.768	0	0th order correction for pk
PHC1 [degrees]	-46.451	0	1st order correction for pk
PH_mod	pk	mc	Phasing modes for trf, xfb, ...

- 13) Type **xf2** (return) to transform the 2D data set followed by **abs2** (return) to baseline smooth it.
- 14) Provided phase constants have been correctly set up the transformed data set plot should resemble that shown below. Black = a positively phased signal, red = a negatively phased signal.




3.0 T₁ Data Set Processing


- 15) Open Topspin's **Analysis** menu and click its **T1/T2 Relaxation** tab. If other Bruker processing software has been installed on the spectrometer terminal, as may be the case on the AVII-600, select the **Analysis** menu's **Topspin T1/T2 module** tab and open its **T1/T2 Relaxation** sub-menu tab.




This will open up the screen display shown below.




- 16) Click the **Extract Slice** button and then click the **Spectrum** button in the panel that appears.




- 17) Enter **1** in the **Slice Number** cell and click the **OK** button. This will display the transformed spectrum ex the first (longest) **VDLIST** value = the one that was phased via the **rser 1** routine in steps 8 and 9.




- 18) Expand the spectrum that appears using Topspin's "E" (Expand) menu button and integrate selected peaks in it as per standard ^1H spectrum processing. Integrals should start and terminate as close as possible to the edge of a peak.




- 19) Click the "Save Region As" icon button (= 3rd from the right in TS3.2's the upper menu bar, or towards the left hand side of the menu bar in TS2.1) = *the one with the floppy disk icon + A below it*) and then its **Export Regions to Relaxation Module and .ret. tab.**





- 20) Click the **Relaxation Window** button, note any screen messages that may appear and close them. Select **Intensity** or **Area** in the plot window that appears. **Area** is often the better choice.



- 21) Click the **Fitting Function** and note the comments about **intensity** or **area** options in the screen display that appears and close it.




- 22) When the Fitting Function message screen is closed a panel with relaxation parameters (split into two screen captures below) will appear.



- 23) Check **Function Type** = **uxnmrt1** and **List file name** = **vdlist**.
Other cells/values can be left as they are (= default settings)

.... next page

- 24) Click the **double red arrow** in the menu bar at the top left of the plot window and **NOT** the single red arrow button in the NMR Relaxation Guide menu below the Fitting Function button.



- 25) The + and - buttons in the upper menu bar can be used to advance (or reverse) the individual T₁ plots.
- 26) Fitting type (**Intensity** or **Area**) can be changed in the plot display. If this is done clicking the **double red arrow** will recalculate the T₁ results and update their plots.
- 27) The NMR Relaxation Guide has buttons which can be used to display and/or print T₁ results. A sample report for one signal is shown below.

1	Dataset :
2	C:\Bruker\TopSpin3.0\data\AVII-600 Set up\nmr\AVII T1+T2/3/pdata/1
3	AREA fit :
4	$I[t] = I[0] + P \cdot \exp(-t/T1)$
5	
6	8 points for Integral 1, Integral Region from 8.708 to 8.668 ppm
7	Results Comp. 1
8	
9	I[0] = 1.093e+000
10	P = -2.024e+000
11	T1 = 3.349s
12	SD = 1.021e-002
13	
14	tau ppm integral intensity
15	
16	10.000s 8.678 5.5583e+008 6.5345e+007
17	6.000s 8.678 4.1279e+008 4.7869e+007
18	3.000s 8.678 1.4514e+008 1.4979e+007
19	1.000s 8.678 -2.1925e+008 -3.0296e+007
20	500.000m 8.678 -3.5543e+008 -4.7347e+007
21	250.000m 8.678 -4.3437e+008 -5.7023e+007
22	100.000m 8.678 -4.8692e+008 -6.3715e+007
23	30.000m 8.678 -5.153e+008 -6.7561e+007