



KJM 9250

Carbon T₁ Spectra on the AVI-600 and AVII-600

Version 1.0



© Professor Emeritus Alistair Lawrence Wilkins,
University of Waikato, New Zealand.
February 2018

Carbon T₁ Spectra on the AVI-600 and AVII-600

1.0 Introduction

An **awcarbont1** parameter set and linked **VDLIST** file have been set up on the AVI-600 and AVII-600 spectrometers running under TS2.1 and TS3.2 respectively.

The **awcarbont1** parameter set has **D1 = 10 sec** and a linked **VDLIST** file with 8 x delays in the range 0.03 to 10 sec with the longest delay (10 sec) the first used VDLIST value and the shortest delay (0.03 sec) the last used value.

D1 and the longest (first) **VDLIST** value should be 3-5 times the longest T₁ in the sample compound.

The **awcarbon16t1** parameter set with a linked **VDLIST** file with **VD_{max}** and **D1 = 16 sec** can be used for signals with longer T₁'s.


T₁ data sets can be processed on the spectrometer terminal or off line using any version of Topspin provided a copy of the linked **VDLIST** file is resident in the experiment's top level folder (= default set up: see below) or a copy of the originally named **VDLIST** file is recreated in the offline terminal's *C:\Bruker\Topspin...\exp\stan\nmr\lists\vd* folder.

 uxnmr.info	2/01/2017 10:50 p...	INFO File	4 KB
 uxnmr.par	2/01/2017 10:50 p...	PAR File	24 KB
 vdlst	10/01/2018 8:49 a...	File	1 KB


NB: The experimental copy of the variable delay file is named as **vdlst** irrespective of the name of the source vdlst file.

2.0 Experiment Set Up


- 1) Create an experiment and read in the **awcarbont1** parameter set (+ **getprosol**).
Pulse programme = **t1irpg** (AVI-600) or **awt1irpg** (AVII-600).
- 2) Review default settings. These settings can be adjusted if required.
TD(F2) = 32K or 64K, TD(F1) = 8.
SW = 240 ppm, O1P = 110 ppm.
D1 = 10 sec.
NS = multiple of 4 or 8, DS = 4 or 8.
- 3) Type **ased** (enter) and review other parameters used in the job including the linked **VDLIST** file = **AWCARBONT1-8**
- 4) The **VDLIST** file should have the entries shown below (next page). Values are in seconds. Do not alter the values in this file. An alternatively named **VDLIST** file should be created if different **VD** and **D1** values are required for a particular compound.



- 5) Set receiver gain using **RGA** (*important!*).
- 6) Type **edp** (enter) and check that **SI(F2) = 32K or 64K**, **SI(F1) = 8**, **WDW(F2) = EM**, **LB (F2) = 2 Hz** or other value of your choice.



- 7) Start the acquisition using the **ZG** command.
- 8) When the experiment has run type **rser 1** (enter) to read in the first serial file which will appear in a TEMP screen display window.
- 9) Type **EFP** (return) to transform it and phase it as per a normal carbon spectrum. The **multiabsn** command with **n = 30-40** can be used to straighten the baseline.



- 10) Type **edp** (enter) and note the phase constants for this spectrum. AVI and AVII-600 cyroprobe ^{13}C spectra will have a large negative **PHC1** value.

PHC0 [degrees]	171.445
PHC1 [degrees]	-856.632
PH_mod	pk


- 11) Close the **TEMP** window and reload the T_1 data set file.

- 12) Type **edp** (enter) and enter the phase constants noted in step 10 above into the **F2 PHC0** and **F2 PHC1** cells and check **PH_MOD = pk**. **F1** cell info is not used.

Phase correction			
PHC0 [degrees]	171.445	0	0th order correction for
PHC1 [degrees]	-856.632	0	1st order correction for
PH_mod	pk	mc	Phasing modes for trf, x


- 13) Type **xf2** (return) to transform the 2D data set followed by **abs2** (return) to baseline smooth it.

- 14) Provided phase constants have been correctly set up the transformed data set plot should resemble that shown below. Black = a positively phased signal, red = a negatively phased signal.




3.0 T₁ Data Set Processing


- 15) Open Topspin's **Analysis** menu and click its **T1/T2 Relaxation** tab. If other Bruker processing software has been installed on the spectrometer terminal, as may be the case on the AVII-600, select the **Analysis** menu's **Topspin T1/T2 module** tab and open its **T1/T2 Relaxation** sub-menu tab.




This will open up the screen display shown below.




- 16) Click the **Extract Slice** button and then click the **Spectrum** button in the panel that appears.




- 17) Enter **1** in the **Slice number** cell and click the OK button. This will display the transformed spectrum ex the first (longest) **VDLIST** value = the one that was phased via the **rser 1** routine in steps 7 and 8.




- 18) Expand the spectrum that appears using Topspin's "E" (Expand) menu button and integrate selected peaks in it as per standard ¹H spectrum processing. Integrals should start and terminate as close as possible to the edge of a peak.




- 19) Click the "**Save Region As**" icon button (= 3rd from the right in TS3.2's the upper menu bar, or towards the left hand side of the menu bar in TS2.1) = *the one with the floppy disk icon + A below it*) and then its **Export Regions to Relaxation Module and .ret. tab**.




- 20) Click the **Relaxation Window** button, note any screen messages that may appear and close them. Select **Intensity** or **Area** in the plot window that appears. **Area** is often the better choice.



- 21) Click the **Fitting Function** button and note the comments about **intensity** or **area** options in the screen display that appears and close it.




- 22) When the Fitting Function message screen is closed a panel with relaxation parameters (split into two screen captures below) will appear.



- 23) Check **Function Type = uxnmrt1** and **List file name = vdlist**. Other cells/values can be left as they are (= default settings)

- 24) Click the **double red arrow** in the menu bar at the *top left* of the plot window and **NOT** the single red arrow button in the NMR Relaxation Guide menu below the Fitting Function button.



- 25) The + and - buttons in the upper menu bar can be used to advance (or reverse) the individual T₁ plots.
- 26) Fitting type (**Intensity** or **Area**) can be changed in the plot display. If this is done clicking the **double red arrow** will recalculate the T₁ results and update their plots.
- 27) The NMR Relaxation Guide has buttons which can be used to display and/or print T₁ results. A sample report for one signal is shown below.

```

1 Dataset :
2 C:\Bruker\TopSpin3.0\data\AVII-600 Set up\nmr\AVII T1+T2/5/pdata/1
3 AREA fit :
4 I[t]=I[0]+P*exp(-t/T1)
5
6 8 points for Integral 1, Integral Region from 131.850 to 131.557 ppm
7 Results      Comp. 1
8
9 I[0] = 9.749e-001
10 P = -1.761e+000
11 T1 = 386.763m
12 SD = 2.981e-002
13
14 tau ppm integral intensity
15
16 10.000s 131.600 1.7524e+008 7.4892e+007
17 6.000s 131.600 1.7987e+008 7.4859e+007
18 3.000s 131.600 1.813e+008 7.4184e+007
19 1.000s 131.600 1.4143e+008 6.6455e+007
20 500.000m 131.600 9.2506e+007 3.7304e+007
21 250.000m 131.601 1.5134e+007 eliminated
22 100.000m 131.600 -7.5088e+007 -2.8382e+007
23 30.000m 131.600 -1.1773e+008 -5.1553e+007

```