

DPX 300 Selective 1D experiments.
(SELNOESY and SELROESY does not work on DPX 300.)
SELCOSY

1. edc define name user expno and so forth
2. rpar proton all
3. getprosol
4. rga
5. zg
6. efp, apk, abs
7. Note down the O1 value of peak to be irradiated
8. edc increase expon with 1
9. rpar uioselcosy all (do NOT type getprosol)
10. Enter the O1 value of the peak to irradiate
11. rga
12. ns and TD0 to be set (total number of scans = NS x TD0) usually NS = 16.
13. zg
14. The experiment can be stopped by typing STOP (not HALT)
15. Do not type tr while doing the experiment
16. efp (NB efp antiphase correlation peaks)
17. Alternatively: FT, abs, PS, /8 several times (all peaks positive).
- 18.

SELTOCSY (better quality than selcosy)

1. edc define name user expno and so forth
2. rpar proton all
3. getprosol
4. rga
5. zg
6. efp, apk, abs
7. Note down the O1 value of peak to be irradiated
8. edc increase expno with 1
9. rpar uioseltocsy5 all (short range correlations) (do NOT type getprosol)
10. or rpar uioseltocsy80 all (medium range correlations)(do NOT type getprosol)
11. or rpar uioseltocsy240 all (long range correlations) (do NOT type getprosol) *This the recommended one!*
12. Enter the O1 value of the peak to irradiate
13. ns and TD0 to be set (total number of scans = NS x TD0) usually NS = 16
14. rga
15. zg
16. The experiment can be stopped by typing STOP (not HALT)
17. Do not type tr while doing the experiment
18. For the 240 and 80 versions: efp (NB all peaks should be positive and phased as for a standard 1H spectrum)
19. For the 5 version
Alternatively: If some peaks show anti phase distortion (with shorter mixing times) power mode processing can be used to generate positive peaks: type: FT, PS, abs, /8 several times (all peaks positive)
19. To get an NMR spectrum with much higher intensity do as follows:
20. For resolution enhancement set LB = -1.5 (or -2), GB = 0.33
21. GFP
22. PS
23. abs
24. /8 several times (all peaks are positive with resolution enhancement processing).
25. LB to 0.1 (or other normal value)
26. GB = 0 BEFORE using EFP to reprocess the spectrum as a conventional 1H spectrum (otherwise EFP with GFP parameters gives rubbish/noise).